

PROCESSING & EMBEDDING ACCESSORIES



Raymond A Lamb's MicroFlow, Slotted and Processing Cassettes have been designed to provide the best performance and reliability when used in our MicroWriter cassette labelers (as sold in the U.S. by Shandon and TBS). These high quality cassettes are available in a variety of colors.

MicroWriter Processing Cassette

MicroWriter Processing Cassettes are available in white, blue, green, pink, yellow, gray, orange and red and are supplied in bags of 1000. White plastic lids are sold separately in bags of 1000 (Cat Ref E10/LID).



Cat Ref	Pack Size	Cat Ref	Pack Size
E10/WHITE	Bag/1000	E10/ORANGE	Bag/1000
E10/BLUE	Bag/1000	E10/RED	Bag/1000
E10/GREEN	Bag/1000	E10/PINK	Bag/1000
E10/YELLOW	Bag/1000	E10/GRAY	Bag/1000
E10/LID	Bag/1000		



MicroWriter Microflow Cassettes

The MicroWriter MicroFlow Cassette is designed to allow better infiltration during processing and to fit most popular base molds for simple embedding. The MicroFlow cassette is available in white, blue, green, pink, yellow, gray, orange and red and supplied in bags of 1000. White plastic MicroFlow lids are sold separately in bags of 1000 (Cat Ref E10/B/LID).

Cat Ref	Pack Size	Cat Ref	Pack Size	Cat Ref	Pack Size
E10/B/WHITE	Bag/1000	E10/B/BLUE	Bag/1000	E10/B/GREEN	Bag/1000
E10/B/PINK	Bag/1000	E10/B/YELLOW	Bag/1000	E10/B/GRAY	Bag/1000
E10/B/ORANGE	Bag/1000	E10/B/RED	Bag/1000	E10/B/LID	Bag/1000

MicroWriter Slotted Cassettes

The MicroWriter Slotted Cassettes are available in white, blue, green, pink, yellow, gray, orange and red and are supplied in bags of 1000. White plastic slotted lids are sold separately in bags of 1000 (Cat Ref E10/S/LID).

Cat Ref	Pack Size	Cat Ref	Pack Size	Cat Ref	Pack Size
E10/S/WHITE	Bag/1000	E10/S/BLUE	Bag/1000	E10/S/GREEN	Bag/1000
E10/S/PINK	Bag/1000	E10/S/YELLOW	Bag/1000	E10/S/GRAY	Bag/1000
E10/S/ORANGE	Bag/1000	E10/S/RED	Bag/1000	E10/S/LID	Bag/1000

METAL CASSETTE LIDS

Permanently reusable perforated stainless steel self-locking covers for processing cassettes.

Cat Ref	Pack Size
E10.85	Pack/25



CASSETTE DIVIDER

A stainless steel cross that fits into a standard processing or embedding cassette to divide the cassette chamber into four sections. Sold individually.

Cat Ref
E10.86



CELLSAFE™

The CellSafe is a disposable plastic biopsy insert consisting of two hinged interlocking frames with an integral nylon mesh. The mesh is extra fine and therefore ensures a high level of specimen security.

Cat Ref	Pack Size
E10.87	Pack/100





BIOPSY INSERT

Two part loose fitting stainless steel capsule with welded mesh which fits inside processing cassettes preventing the loss of small or friable material while allowing interchange of fluids. Sold individually.

Cat Ref	Internal Size	Mesh Hole Size
E10.95	23mm diameter x 3.5mm deep	0.5mm



CURETTING CASSETTES

Modified MicroWriter processing cassette with mesh welded to a steel lid and mesh insert. Prevents the loss of small specimens while allowing fluid interchange. Sold individually.

Cat Ref	Internal Size	Mesh Hole Size
E10.9	30 x 24 x 3.5mm deep	0.5mm



BIOPSY INSERT PADS

Foam pad inserts cut to the size of a processing cassette chamber available in blue or gray. Biopsy material is placed between two pads in the cassette and can be processed in the normal way.

Cat Ref	Size	Pack Size
E10.75/BLUE	26 x 31 x 2mm deep	Bag/1000
E10.75/GRAY	26 x 31 x 2mm deep	Bag/1000



CASSETTE INSERT PAPER

Whatman qualitative cellulose paper cut to the size of a processing cassette chamber. Biopsy material is placed between two sheets in the cassette and can be processed in the normal way.

Cat Ref	Size	Pack Size
E10.7	26 x 31mm	Pack/1000



METAL BASE MOLDS

Single reusable stainless steel molds for processing cassettes.

Cat Ref	Mold Size	Pack Size
E10/M1	7 x 7 x 5mm deep	Box/5
E10/M2	24 x 24 x 5mm deep	Box/5
E10/M3	30 x 24 x 5mm deep	Box/5
E10/M4	37 x 24 x 5mm deep	Box/5



DISPOMOLDS FOR PROCESSING CASSETTES

Disposable plastic molds used with processing cassettes. Each strip contains 5 molds.

Cat Ref	Mold Size	Pack Size
E10.8/L7	7 x 7 x 5mm deep	100 strips
E10.8/L15	15 x 15 x 5mm deep	100 strips
E10.8/L24	24 x 24 x 5mm deep	100 strips
E10.8/L30	30 x 24 x 5mm deep	100 strips
E10.8/L37	37 x 24 x 5mm deep	100 strips



DISPOMOLDS FOR EMBEDDING RINGS

Disposable plastic molds used with embedding rings. Each strip contains 5 molds.

Cat Ref	Mold Size	Pack Size
E10.6/2001/33	33 x 25 x 13mm deep	100 strips
E10.6/2001/39	39 x 25 x 13mm deep	100 strips



EMBEDDING RINGS

Sturdy plastic molding that forms the permanent base of a paraffin block when used with Dispomolds. Provides a 10mm deep base to clamp to a microtome.

Cat Ref	Dimensions	Pack Size
E10.6	41 x 28 x 10mm deep	Pack/500