

LAMB'S LABSTACK FILING CABINETS

Raymond A Lamb manufactures a wide variety of high-quality, epoxy-coated, stackable steel cabinets for the storage of microscope slides, slide filing wallets, processing cassettes, tissue blocks, transparencies and A5 index cards.

The range is available in two sizes and can be adapted to meet any storage requirements.

The standard size for floor stacking is available in three colors: BLUE, GREEN and TWO-TONE GRAY.

Each cabinet is completely sealed top and bottom to ensure dust-free storage and has a lip running the full width of the front and rear to allow it to fit securely with the next cabinet in the stack. Each drawer is fitted with stops to prevent accidental removal, however, the drawers can be removed and relocated.

Base stands (black only) and colored top covers are available for the cabinets. The base stand comes in a standard model and a heavy-duty wheeled base model and the top covers allow cabinets to also become additional work surfaces.

The following pages provide full details of each cabinet and the accessories available.

Non standard cabinets can be manufactured to order – please contact Raymond A Lamb for details.





BLOCKSTORE CABINET

The Blockstore cabinet has eight drawers and accepts embedding rings, processing cassettes and unmounted blocks up to 12 x 25 x 30mm. Each drawer has dividers producing seven rows 30mm wide. Unmounted blocks are placed in cardboard boxes as illustrated.

| Cat Ref | Capacity Processing Cassettes | Capacity E43.23 | Capacity Embedding Rings | Dimensions | Weight (empty) |
|--------------------------|-------------------------------|-----------------|--------------------------|-------------------|----------------|
| E43.22/BLUE, GREEN, GRAY | 3000 | 2050 | 1230 | 480 x 480 x 275mm | 17.5kg |
| E43.22/MINI (Gray only) | 1400 | 1000 | 600 | 480 x 240 x 275mm | 11Kg |

NB Capacities for embedding rings and processing cassettes are for typical blocks not the mounts alone.

ACCESSORIES FOR BLOCKSTORE CABINETS

CARDBOARD BOXES FOR UNMOUNTED BLOCKS


| Cat Ref | Dimensions | Pack |
|---------|----------------|------|
| E43.23 | 13 x 26 x 32mm | 1000 |



CASSETTE STORAGE BOX

Stainless steel box holds approximately 200 processing cassettes in 5 rows.

| Cat Ref | Dimensions |
|----------|------------------|
| E43.22/B | 310 x 162 x 46mm |



7 - DRAWER FILING CABINET

The 7 - Drawer Filing Cabinet can be used either for the storage of larger microscope slides up to 76 x 51mm or paraffin wax blocks etc in cardboard boxes and Super Mega Cassettes. It may also be used for the storage of photographic transparencies and most grid boxes used in electron microscopy. If a drawer is fitted with our E43.21 divider two rows of 76 x 26 slides can be stored.

| Cat Ref | Capacity Slides | Capacity E43.28 | Capacity E43.29 | Dimensions | Weight (empty) |
|--------------------------|-----------------|-----------------|-----------------|-------------------|----------------|
| E43.13/BLUE, GREEN, GRAY | 3000 | 133 | 49 | 480 x 480 x 140mm | 10Kg |
| E43.13/MINI (Gray only) | 1500 | 63 | 21 | 480 x 240 x 140mm | 5.5Kg |

ACCESSORIES FOR 7-DRAWER CABINET

CARDBOARD BOXES FOR UNMOUNTED BLOCKS

| Cat Ref | Dimensions | Pack |
|---------|------------------|------|
| E43.28 | 49.5 x 25 x 80mm | 100 |
| E43.29 | 63 x 50 x 80mm | 100 |

DRAWER DIVIDER

Stainless steel divider fits into 7-Drawer Cabinet drawer dividing it to allow for storage drying.

| Cat Ref | For Cabinet |
|-------------|-------------|
| E43.21 | E43.13 |
| E43.21/MINI | E43.13/MINI |

SPRING SPACER

Fits into 7-Drawer Cabinet Drawer to separate transparencies and slides while drying.

| Cat Ref | For Cabinet |
|------------|-------------|
| E43.17 | E43.13 |
| E43.17MINI | E43.13/MINI |

14 - DRAWER SLIDE FILING CABINET

The 14-Drawer Cabinet can be used for the storage of mainly 76 x 26mm microscope slides vertically.

| Cat Ref | Capacity | Dimensions | Weight (empty) |
|--------------------------|-------------|-------------------|----------------|
| E43.12/BLUE, GREEN, GRAY | 6000 slides | 480 x 480 x 140mm | 11kg |
| E43.12/MINI (Gray only) | 3000 slides | 480 x 240 x 140mm | 6.7Kg |



ACCESSORIES FOR 14-DRAWER CABINET

SPRING SPACER

Fits into 14-Drawer Cabinet drawer to separate slides while drying.

| Cat Ref | For cabinet |
|-------------|-------------|
| E43.16 | E43.12 |
| E43.16/MINI | E43.12/MINI |



SLIDE FILING LOCATION SLIPS

Detailed colored slips for the replacement of borrowed slides. Available in booklets of 100 slips.

| Cat Ref | Color |
|---------|-------|
| E43.25 | White |
| E43.26 | Pink |
| E43.27 | Green |



HORIZONTAL SLIDE FILING CABINET

The Horizontal cabinet has two drawers with polypropylene dividers grooved to hold 76 x 26mm microscope slides horizontally.


| Cat Ref | Capacity | Dimensions | Weight (empty) |
|--------------------------|-------------|-------------------|----------------|
| E43.33/BLUE, GREEN, GRAY | 1230 slides | 480 x 480 x 140mm | 11kg |
| E43.33/MINI (Gray only) | 600 slides | 480 x 240 x 140mm | 6.5Kg |



STEEL SLIDE BOX

Stainless steel box with fitted lid and card index holder. It has internal dividers for 76 x 26mm microscope slides.

| Cat Ref | Capacity | Dimensions |
|---------|------------------|------------------|
| E43.11 | Over 1000 slides | 312 x 153 x 83mm |





SLIDE WALLET FILING CABINET

The Slide Wallet Cabinet has three drawers for the storage of Slide Wallets (Cat Ref E48.5).

| Cat Ref | Capacity | Dimensions | Weight (empty) |
|--------------------------|-------------|-------------------|----------------|
| E43.10/BLUE, GREEN, GRAY | 1800 slides | 480 x 480 x 140mm | 8kg |
| E43.10/MINI (Gray only) | 900 slides | 480 x 240 x 140mm | 5.5Kg |



SLIDE FILING WALLET FOR E43.10 CABINET

Cardboard wallet to take 76 x 26mm microscope slides.


| Cat Ref | Capacity | Size | Pack |
|---------|----------|------------|------|
| E48.5 | 3 slides | 85 x 126mm | 100 |



INDEX CARD CABINET

The Index Card Cabinet has two drawers for the storage of index cards (holds A5 and 5" x 8" cards). The capacity of each drawer depends on the type of card used. Available in GRAY only.

| Cat Ref | Dimensions | Weight (empty) |
|----------------------------|-------------------|----------------|
| E43.10/A5/GRAY | 480 x 480 x 140mm | 8 Kg |
| E43.10/A5/MINI (Gray only) | 480 x 240 x 140mm | 5.5Kg |



TOP COVER

Top covers are available in BLUE, GREEN and TWO-TONE GRAY for the standard size cabinets. A MINI Top Cover is in a matching GRAY.

| Cat Ref | Dimensions | For Cabinet |
|--------------|------------------|----------------------------|
| E43.14/BLUE | 480 x 480 x 25mm | All standard size cabinets |
| E43.14/GREEN | 480 x 480 x 25mm | All standard size cabinets |
| E43.14/GRAY | 480 x 480 x 25mm | All standard size cabinets |
| E43.14/MINI | 480 x 240 x 25mm | All MINI cabinets |



BASE STAND

Black epoxy coated base stand.

| Cat Ref | Dimensions | For Cabinet |
|-------------|-------------------|----------------------------|
| E43.15 | 480 x 480 x 110mm | All standard size cabinets |
| E43.15/MINI | 480 x 240 x 110mm | All MINI cabinets |



HEAVY DUTY WHEEL BASE STAND

Black epoxy coated reinforced base stand incorporating four heavy-duty swivel castor wheels. Only available for the standard size cabinets.

| Cat Ref | Dimensions | Cabinets |
|-------------|-------------------|----------------------------|
| E43.15/HDWB | 480 x 240 x 200mm | All standard size cabinets |

SLIDE MAILING BOXES

Heavy-duty cardboard box with cloth covered edges and a fitted lid. The wooden slide rack inside holds 76 x 26mm microscope slides.

| Cat Ref | Capacity |
|---------|-----------|
| E33 | 2 slides |
| E34 | 4 slides |
| E35 | 6 slides |
| E36 | 12 slides |
| E37 | 25 slides |



HINGED SLIDE BOXES

Leatherette covered card box with metal hinges and clasp. The wooden slide rack inside holds 76 x 26mm microscope slides.


| Cat Ref | Capacity |
|---------|------------|
| E32.5 | 50 slides |
| E32.6 | 100 slides |



PLASTIC SLIDE BOXES

Molded high-impact polystyrene box with a close-fitting snap-on lid and separate card index. It is grooved to take 76 x 26mm microscope slides.

| Cat Ref | Capacity |
|---------|------------|
| E38 | 25 slides |
| E39 | 100 slides |



HINGED PLASTIC SLIDE BOX

Silicone coated black polystyrene box with plastic hinge and catch and an index inside the lid. It is recessed to take 76 x 26mm microscope slides.

| Cat Ref | Capacity |
|---------|----------------------------------|
| E40.5 | 100 slides (200 back to back) |



WOODEN SLIDE BOXES

Polished wooden box with brass hinge and catch and an index inside the lid. It is grooved to take 76 x 26mm microscope slides.

| Cat Ref | Capacity |
|---------|------------|
| E39.5 | 50 slides |
| E40 | 100 slides |





2-SLIDE MAILER

Molded polythene with attached snap-on lid and grooves to hold 76 x 26mm microscope slides.

| Cat Ref | Capacity | Weight* |
|---------|----------|---------|
| E6.2 | 2 slides | 20g |



5-SLIDE MAILER

Molded polypropylene with attached snap-on lid and grooves to hold five 76 x 26mm microscope slides.

| Cat Ref | Capacity | Weight* |
|---------|----------|---------|
| E6.3 | 5 slides | 33g |



SLIDE MAILER TRAY

Heavy duty cardboard construction with cloth hinge. The tray is recessed to take 76 x 26mm microscope slides and may be mailed without additional packaging using labels E48.1 illustrated below.

| Cat Ref | Capacity | Weight* |
|---------|----------|---------|
| E48 | 2 slides | 40g |



MAILING LABELS

Self-adhesive labels printed "FRAGILE-MICROSCOPE SLIDES". For use with E48 Slide Mailer Trays or in general packing.

| Cat Ref | Pack |
|---------|------|
| E48.1 | 50 |

* Weights are calculated to include the full capacity of each container using 1.0-1.2mm thickness microscope slides.

SLIDE TRAY BOX

The box opens at the front and has a fitted lid.

| Cat Ref | Dimensions | Capacity |
|---------|------------------|----------|
| E32 | 265 x 112 x 68mm | 12 Trays |

9 - SLIDE TRAY

Heavy duty cardboard construction with cloth hinged flaps and a single compartment holding slides of various sizes.

| Cat Ref | Compartment Size | Capacity |
|---------|------------------|----------|
| E49 | 235 x 78mm | 9 slides |



20 - SLIDE TRAY

Heavy duty cardboard tray with cloth hinged flaps completely covering the slides. The tray is recessed to accept 76 x 26mm slides individually with cut-away for easy removal.

| Cat Ref | Dimensions | Capacity |
|---------|-------------|-----------|
| E44 | 200 x 350mm | 20 slides |



24 - SLIDE TRAY

Heavy duty cardboard tray with cloth hinged flaps completely covering the slides. The tray has four compartments to accept microscope slides of various sizes. Compartment size 82 x 156mm.

| Cat Ref | Dimensions | Capacity |
|---------|-------------|-----------|
| E45 | 200 x 350mm | 24 slides |



SLIDE TRAY RACKS & ALUMINUM SLIDE TRAYS

Stainless steel slide tray racks to accept aluminum slide trays. The rack allows air to circulate over the slides and both sizes of racks hold up to 100 slides. The aluminum slide trays have elevated ridges to separate 76 x 26mm slides and finger holes for easy removal. Ideal for drying fresh mounts.

| Cat Ref | | Dimensions | Capacity |
|----------|-------------------------------------|-------------------|-------------------|
| E45.SR10 | Slide Tray Rack | 330 x 212 x 112mm | 10 x E45/10 trays |
| E45.SR5 | Slide Tray Rack | 330 x 197 x 112mm | 5 x E45/20 trays |
| E45/10 | Aluminum Slide Tray (fits E45/SR10) | 330 x 100mm | 10 slides |
| E45/20 | Aluminum Slide Tray (fits E45/SR5) | 330 x 190mm | 20 slides |

