



### FREEZER INVENTORY RANGE

The freezer inventory range is manufactured from polished stainless steel and is designed to get the maximum storage capacity using standard units. Units with different shelf heights can also be incorporated. The design will also enable the user to reorganize or add to existing freezer storage without having to replace the entire system.



### FRONT OPENING FREEZER RANGE

This range is designed to stack and lock together. Shelf units are available with 2, 3, 4, 5 or 6 shelf spaces in a range of standard sizes. The shelves are made up of rows of 6mm tabs pushed in from the sides.

The picture above shows how the units stack and lock together using a small 2mm tab at the top of each unit and a corresponding notch at the bottom.



### CHEST FREEZER RANGE

Shelf units for this range are available with 5, 6, 7, 8, 9 or 10 shelf spaces in a range of standard sizes. The shelves are made up of rows of 6mm tabs pushed in from the sides. This makes the units much more cost effective.

### MICROTITRATION PLATE UNIT

The Microtitration plate units are available in both front opening and chest freezer configurations. Front opening units are supplied with stainless steel trays – the chest freezer version does not require trays.



### TRAYS FOR FRONT OPENING AND CHEST SHELF UNITS

Trays can be supplied empty, with a range of colored divisions or with an aluminum rack as shown above. Trays are ordered separately to the shelf units to make it easy to mix storage options.

### NON-STANDARD SHELF UNITS AND TRAYS - "SPECIALS"

Using our new manufacturing process we are happy to supply non-standard shelf units and trays to suit individual requirements. Prices for non-standard or "special" units will be based on our standard range.

# FRONT OPENING FREEZER SHELF UNITS



To order shelf units you need to know the available width, depth and height of your freezer. Then select units to make up the required width with a depth suitable for your freezer from the first column. You can then decide the combination of shelves and units from the central column to best utilize the height available in your freezer. To order simply quote the catalog reference in the right hand column.

Units with Shelf Spacing 56mm			Cat Ref
140mm Wide x 385mm Deep	2 shelves	118 mm High	EU140W385D56S2
	3 shelves	176 mm High	EU140W385D56S3
	4 shelves	234 mm High	EU140W385D56S4
	5 shelves	292 mm High	EU140W385D56S5
	6 shelves	350 mm High	EU140W385D56S6
	160mm Wide x 385mm Deep	2 shelves	118 mm High
3 shelves		176 mm High	EU160W385D56S3
4 shelves		234 mm High	EU160W385D56S4
5 shelves		292 mm High	EU160W385D56S5
6 shelves		350 mm High	EU160W385D56S6
140mm Wide x 460mm Deep		2 shelves	118 mm High
	3 shelves	176 mm High	EU140W460D56S3
	4 shelves	234 mm High	EU140W460D56S4
	5 shelves	292 mm High	EU140W460D56S5
	6 shelves	350 mm High	EU140W460D56S6
	160mm Wide x 460mm Deep	2 shelves	118 mm High
3 shelves		176 mm High	EU160W460D56S3
4 shelves		234 mm High	EU160W460D56S4
5 shelves		292 mm High	EU160W460D56S5
6 shelves		350 mm High	EU160W460D56S6
140mm Wide x 550mm Deep		2 shelves	118 mm High
	3 shelves	176 mm High	EU140W550D56S3
	4 shelves	234 mm High	EU140W550D56S4
	5 shelves	292 mm High	EU140W550D56S5
	6 shelves	350 mm High	EU140W550D56S6
	160mm Wide x 550mm Deep	2 shelves	118 mm High
3 shelves		176 mm High	EU160W550D56S3
4 shelves		234 mm High	EU160W550D56S4
5 shelves		292 mm High	EU160W550D56S5
6 shelves		350 mm High	EU160W550D56S6
Units with Shelf Spacing 70mm			Cat Ref
140mm Wide x 385mm Deep	2 shelves	146 mm High	EU140W385D70S2
	3 shelves	218 mm High	EU140W385D70S3
	4 shelves	290 mm High	EU140W385D70S4
	5 shelves	362 mm High	EU140W385D70S5
	160mm Wide x 385mm Deep	2 shelves	146 mm High
3 shelves		218 mm High	EU160W385D70S3
4 shelves		290 mm High	EU160W385D70S4
5 shelves		362 mm High	EU160W385D70S5
140mm Wide x 460mm Deep		2 shelves	146 mm High
	3 shelves	218 mm High	EU140W460D70S3
	4 shelves	290 mm High	EU140W460D70S4
	5 shelves	362 mm High	EU140W460D70S5
	160mm Wide x 460mm Deep	2 shelves	146 mm High
3 shelves		218 mm High	EU160W460D70S3
4 shelves		290 mm High	EU160W460D70S4
5 shelves		362 mm High	EU160W460D70S5
140mm Wide x 550mm Deep		2 shelves	146 mm High
	3 shelves	218 mm High	EU140W550D70S3
	4 shelves	290 mm High	EU140W550D70S4
	5 shelves	362 mm High	EU140W550D70S5
	160mm Wide x 550mm Deep	2 shelves	146 mm High
3 shelves		218 mm High	EU160W550D70S3
4 shelves		290 mm High	EU160W550D70S4
5 shelves		362 mm High	EU160W550D70S5

Units with Shelf Spacing 56mm			Cat Ref	
139mm Wide x 143mm Deep	5 shelves	301 mm High	EC139W143D56S5	
	6 shelves	359 mm High	EC139W143D56S6	
	7 shelves	417 mm High	EC139W143D56S7	
	8 shelves	475 mm High	EC139W143D56S8	
	9 shelves	533 mm High	EC139W143D56S9	
	10 shelves	591 mm High	EC139W143D56S10	
	139mm Wide x 163mm Deep	5 shelves	301 mm High	EC139W163D56S5
		6 shelves	359 mm High	EC139W163D56S6
		7 shelves	417 mm High	EC139W163D56S7
		8 shelves	475 mm High	EC139W163D56S8
9 shelves		533 mm High	EC139W163D56S9	
159mm Wide x 163mm Deep	5 shelves	301 mm High	EC159W163D56S5	
	6 shelves	359 mm High	EC159W163D56S6	
	7 shelves	417 mm High	EC159W163D56S7	
	8 shelves	475 mm High	EC159W163D56S8	
	9 shelves	533 mm High	EC159W163D56S9	
10 shelves	591 mm High	EC159W163D56S10		

Units with Shelf Spacing 70mm			Cat Ref
139mm Wide x 143mm Deep	5 shelves	371 mm High	EC139W143D70S5
	6 shelves	443 mm High	EC139W143D70S6
	7 shelves	515 mm High	EC139W143D70S7
	8 shelves	587 mm High	EC139W143D70S8
	9 shelves	659 mm High	EC139W143D70S9
139mm Wide x 163mm Deep	5 shelves	371 mm High	EC139W163D70S5
	6 shelves	443 mm High	EC139W163D70S6
	7 shelves	515 mm High	EC139W163D70S7
	8 shelves	587 mm High	EC139W163D70S8
	9 shelves	659 mm High	EC139W163D70S9
159mm Wide x 163mm Deep	5 shelves	371 mm High	EC159W163D70S5
	6 shelves	443 mm High	EC159W163D70S6
	7 shelves	515 mm High	EC159W163D70S7
	8 shelves	587 mm High	EC159W163D70S8
	9 shelves	659 mm High	EC159W163D70S9

Trays are available empty or with divisions made up of interlocking combs (see Specimen Storage page 51 for more details)

Division Sizes  
 10.5mm    25mm  
 12.5mm    29mm  
 16.5mm    32mm  
 18.5mm    35mm  
 22.5mm

Specify division size when ordering

### TRAYS FOR FRONT OPENING FREEZER SHELF UNITS

To Fit Shelf Units	External Dimensions	Cat Ref
140mm Wide x 385mm Deep Units	135mm Wide x 385mm Deep	ET135W385/E or Div size
	135mm Wide x 395mm Deep	ET135W395/E or Div size
140mm Wide x 460mm Deep Units	135mm Wide x 460mm Deep	ET135W460/E or Div size
	135mm Wide x 470mm Deep	ET135W470/E or Div size
140mm Wide x 550mm Deep Units	135mm Wide x 550mm Deep	ET135W550/E or Div size
	135mm Wide x 560mm Deep	ET135W560/E or Div size
160mm Wide x 385mm Deep Units	155mm Wide x 385mm Deep	ET155W385/E or Div size
	155mm Wide x 395mm Deep	ET155W395/E or Div size
160mm Wide x 460mm Deep Units	155mm Wide x 460mm Deep	ET155W460/E or Div size
	155mm Wide x 470mm Deep	ET155W470/E or Div size
160mm Wide x 550mm Deep Units	155mm Wide x 550mm Deep	ET155W550/E or Div size
	155mm Wide x 560mm Deep	ET155W560/E or Div size

### TRAYS FOR CHEST FREEZER SHELF UNITS

To Fit Shelf Units	External Dimensions	Cat Ref
To fit 139mm Wide x 143mm Deep Units	134mm Wide x 134mm Deep	ET134W134/E or Div size
To fit 139mm Wide x 163mm Deep Units	134mm Wide x 154mm Deep	ET134W154/E or Div size
To fit 159mm Wide x 163mm Deep Units	154mm Wide x 154mm Deep	ET154W154/E or Div size