

LAMB SLIDES

High quality, pre-cleaned slides with ground edges, available in plain or twin-frosted varieties. The twin-frosted slides have a finely-etched surface to produce the best marking results. Supplied in boxes of 72.

Cat Ref	Size	Thickness	Finish
E25L	76 x 26mm	1.0mm	Plain
E25.3L	76 x 26mm	1.0mm	Twin Frost



LAMB HISTOBOND SLIDES

Lamb Histobond slides have a permanent, positive charge which secures all types of sections to the slides, thus preventing sections floating off during washing. Sections even adhere to the slides when used in microwave techniques. Supplied in boxes of 100.

Cat Ref	Size	Thickness
E27.5HB	75 x 25mm	1.0mm



COVERSLIPS

Cat Ref E23 Please specify thickness and size when ordering

High quality coverslips, cleaned ready for use. Supplied in ounces.

Circles	Squares	Rectangles	Standard Thicknesses
10mm diameter	18x18mm	22x26mm	No 0 - 0.085-0.130mm
13mm diameter	22x22mm	22x32mm	No 1 - 0.130-0.160mm
16mm diameter	24x24mm	22x40mm	No 1.5 - 0.155-0.185mm
19mm diameter		22x44mm	No 2 - 0.190-0.250mm
22mm diameter		22x50mm	
		22x57mm	Non-Standard Thicknesses
		22x64mm	No 3 - 0.250-0.350mm
		24x32mm	No 5 - 0.470-0.600mm
		24x40mm	
		24x50mm	
		35x64mm	
		48x64mm	
		51x64mm	

OTHER SIZES ARE AVAILABLE - PLEASE CONTACT RAYMOND A LAMB FOR DETAILS



HOG HAIR BRUSHES

Hog hair brushes with metal mount and wooden handle.

Cat Ref	Size
E6	No 4 Medium
E6.1	No 6 Thick



KUEHNE'S COVERSIP TWEEZERS

Stainless steel tweezers with angled flat blades for handling coverslips.

Cat Ref	Length
E16.7	11cm



LENS TISSUE BOOKLETS

Booklet of 50 leaves of Whatman tissues.

Cat Ref	Size
E19	100x75mm



CANADA BALSAM BOTTLES

Glass bottles with loose fitting drop-on cap and glass rod.

Cat Ref	Capacity
E2	30ml
E2.1	50ml



CEDARWOOD OIL BOTTLE

Glass bottle with polypropylene stopper and built-in glass dropping rod.

Cat Ref	Capacity
E5	50ml



POLYSTOP BOTTLES

Glass bottles with plastic stopper and bulb.

Cat Ref	Capacity
E3.1	30ml
E3.2	60ml
E3.3	125ml