



LAMB DPX

The Lamb DPX is a mounting medium ideal for routine use as it is rapid drying. It is supplied in a wide neck, tamper evident UN standard plastic bottle.

Cat Ref LAMB/DPX Pack Size 1/Liter

LAMB LOW VISCOSITY DPX

The Lamb low viscosity formula DPX is ideal for use with automatic coverslipping instruments and is exceptionally fast drying. It is supplied in a wide neck, tamper evident UN standard plastic bottle.

Cat Ref LAMB/DPX/LV Pack Size 1/Liter



LAMB OCT EMBEDDING MEDIUM

The Lamb OCT is an ideal embedding medium for use when cutting tissue sections. It is supplied in a clear plastic bottle with a convenient nozzle molded into the cap.

Cat Ref LAMB/OCT Pack Size 125ml



LAMB IMMERSION OIL

The Lamb Immersion Oil is supplied in a clear plastic bottle.

Cat Ref SCOIL Pack Size 500ml



LAMB COOL-JET

The Lamb Cool-Jet is an environmentally friendly instant freezing aerosol for freezing tissue and paraffin blocks.

Cat Ref C41.5 Pack Size 333g

PAPANICOLAOU STAINS

A range of three stains – EA50 and OG6 counter stains and Harris Hematoxylin – supplied in wide neck, tamper evident UN standard plastic bottles.

- Excellent color balance, clarity and stability
- Each batch is quality controlled by independent laboratories
- Suitable for gynecological as well as non-gynecological applications

	Cat Ref	Pack Size
EA50	LAMB/031-G	4/Quarts
OG6	LAMB/032-G	4/Quarts
Harris Hematoxylin	LAMB/033-G	4/Quarts

Detailed leaflet available on request



HISTOCLEAR / HISTOCLEAR II

Histoclear and Histoclear II are superior non-toxic clearing agents which lead to the production of high quality tissue slides.

	Cat Ref	Pack Size
Histoclear	C78-G	1 Gallon
Histoclear	C78-H	5 Gallon
Histoclear II (reduced citrus odor)	C78/2-G	1 Gallon
Histoclear II (reduced citrus odor)	C78/2-H	5 Gallon

HISTOMOUNT

Histomount is a non-hazardous synthetic mounting medium.

Cat Ref Pack Size

C79-A	100ml
C79-C	450ml



CALCI-CLEAR RAPID

Calci-Clear is a decalcifier ideal for use where speed is important for diagnosis.

Cat Ref Pack Size

C122.5	1 Gallon
--------	----------

WAXES

LambWax Cat Ref W1

Pelletized paraffin wax with added polymers - MP 57- 58°C.
Supplied in packs of 10Kgs (22.5lbs)

PURE PARAFFIN WAX Cat Ref W2

Pelletized paraffin wax specially filtered for highest quality - melting point 56°C.
Supplied in packs of 10Kgs (22.5lbs)

RALWax Cat Ref W11.4

Pelletized pure paraffin wax with plastic co-polymers that give a hard white paraffin of amorphous particle structure - melting point 58 - 60°C.
Supplied in packs of 10Kgs (22.5lbs)

APPLICATIONS

1. Particularly hard or fibrous routine tissue may be cut more easily if, after standard processing in LambWax or Pure Paraffin Wax, they are then vacuum embedded in RALWax.
2. Freeze-dried tissues, where immunohistochemistry is not required, particularly tissues that have been vapor fixed after freeze-drying.
3. Bone sectioning either undecalcified (Sellotape or Tripp and Mackay techniques) or decalcified (e.g. trephine biopsies).

PHYSIOWAX Cat Ref W4

Low melting point wax - melting point 45°C. Supplied in 1Kg (2 x 500g tablets).



DRY DYES / DRY CHEMICALS

	CI No	Cat Ref	Pack Size		CI No	Cat Ref	Pack Size
Acridine Orange	46005	S3-2	25g	Neutral Red	50040	S254-2	25g
Alcian Blue 8GX	74240	S6-2	25g	New Fuchsin	42520	S258-2	25g
Alcian Green 2GX	-	S7-2	25g	New Methylene Blue	52030	S259-2	25g
Alcian Yellow	-	S9-2	25g	Nile Blue (Sulphate)	51180	S263-2	25g
Alizarin Red S	58005	S16-2	25g	Nuclear Fast Red	60760	S265-1	10g
Amido Black	(see Naphthalene Black 10B)			Nuclear Fast Red	60760	S265-2	25g
Aniline Blue (Alcohol Soluble)	42775	S22-2	25g	Oil Red O	26125	S267-2	25g
Aniline Blue (Water Soluble)	42755	S23-2	25g	Orange G	16230	S269-2	25g
Auramine O	41000	S30-2	25g	Orcein Natural	(see S271 Orcein Synthetic as alternative)		
Auramine O	41000	S30-3	100g	Orcein Synthetic	-	S271-0	5g
Azur A (Macneal)	52005	S42-1	10g	Orcein Synthetic	-	S271-2	25g
Biebrich Scarlet	26905	S52-2	25g	Pararosaniline	(see Fuchsin Basic (Feulgen))		
Bismarck Brown Y	21000	S54-2	25g	Phloxine B	45410	S282-2	25g
Brilliant Cresyl Blue	51010	S58-2	25g	Ponceau S	27195	S287-2	25g
Carmine	75470	S70-0	5g	Ponceau de Xylidine	16150	S289-2	25g
Carmine	75470	S70-2	25g	Pontamine Sky Blue	24410	S290-2	25g
Celestin Blue	51050	S75-1	10g	Pyronin Y	45005	S293-1	10g
Celestin Blue	51050	S75-2	25g	Pyronin Y	45005	S293-2	25g
Chrome Cyanine R	43820	S81.5-2	25g	Rhodamine B	45170	S296-2	25g
Chromotrope 2R	16570	S82-2	25g	Saffron	75100	S303-1	10g
Congo Red	22120	S87-2	25g	Safranin O	50240	S305-2	25g
Cresyl Fast Violet	(see S100 Cresyl Violet Acetate as alternative)			Safranin O	50240	S305-3	100g
Cresyl Violet Acetate	-	S100-0	5g	Scarlet R	26105	S306-2	25g
Cresyl Violet Acetate	-	S100-2	25g	Silver Nitrate	-	C128-2	25g
Crystal Ponceau	16250	S101-2	25g	Silver Nitrate	-	C128-3	100g
Crystal Violet	42555	S102-2	25g	Sirius Red F3B	35780	S309-2	25g
Crystal Violet	42555	S102-3	100g	Sudan IV	(see Scarlet R)		
Eosin Ethyl (Alcoholic)	45386	S114-2	25g	Sudan Black B	26150	S318-2	25g
Eosin (Yellowish Water/Alcohol Soluble)	45380	S115-2	25g	Tartrazine	19140	S324-2	25g
Eosin (Yellowish Water/Alcohol Soluble)	45380	S115-3	100g	Thioflavine T	49005	S333-2	25g
Erythrosin B	45430	S122-2	25g	Thionin	52000	S334-2	25g
Fast Green FCF	42053	S142-2	25g	Toluidine Blue	52040	S338-2	25g
Field Stain A	-	S163-2	25g	Victoria Blue B	44045	S345-2	25g
Field Stain B	-	S164-2	25g	Victoria Blue 4R	42563	S346-2	25g
Fuchsin Acid	42685	S168-2	25g	Wright Stain	-	S349-2	25g
Fuchsin Basic (Feulgen)	42500	S169-2	25g				
Fuchsin Basic (Feulgen)	42500	S169-3	100g				
Fuchsin Basic (Ziehl-Neelsen)	42510	S169.5-2	25g				
Fuchsin Basic (Ziehl-Neelsen)	42510	S169.5-3	100g				
Gallocyanine	51030	S170-2	25g				
Giemsa Stain	-	S172-2	25g				
Gold Chloride (Yellow)	-	C76	1g				
Haematein	-	S176-2	25g				
Hematoxylin	75290	S178-2	25g				
Hematoxylin	75290	S178-3	100g				
Hematoxylin (Ehrlich)	-	S180-2	25g				
Hematoxylin (Harris)	-	S181-2	25g				
Jenner Stain	-	S189-2	25g				
Leishman Stain	-	S193-2	25g				
Light Green	42095	S194-2	25g				
Lissamine Fast Yellow	18965	S197-2	25g				
Luxol Fast Blue MBS	-	S205-2	25g				
Malachite Green	42000	S208-2	25g				
Martius Yellow	10315	S210-2	25g				
May Grunwald Stain	-	S212-2	25g				
Methyl Green	42585	S219-2	25g				
Methyl Violet 2B	42535	S223-2	25g				
Methyl Violet 6B	42555	S224-2	25g				
Methylene Blue	52015	S225-2	25g				
Methylene Blue	52015	S225-3	100g				
Methylene Blue (Loeffler Compound)	-	S226-2	25g				
Naphthalene Black 10B	20470	S239-1	10g				
Naphthalene Black 10B	20470	S239-2	25g				