

# Revco -80°C

## Countertop Ultra-Low Temperature Freezer



Model ULT185-5A

Offers  
Full  
1 cu.ft.  
Capacity

### Ideal for Small Storage Requirements

Perfect for clinical trial storage, Revco offers this countertop ultra-low freezer with a full 1 cu.ft. capacity.

### Operating Temperatures

- -80°C, Range, -50°C to -80°C
- Air-cooled refrigeration systems are engineered to maintain design temperatures in warm 32°C (+90°F) ambient conditions

### General Features

- 1 cu. ft. capacity (28 liters)
- Built-in safety alarm system with remote contacts
- Key operated switch for main power
- Digital temperature display resolution to 1°C
- Internal dimensions: 9.63"H x 12" F-B x 15.25"W
- External dimensions: 20"H x 24" F-B x 30.75"W

- Shipping weight 190 lbs

### Cabinet Construction

- Industrial-grade hermetically sealed refrigeration compressor
- Foamed-in-place urethane insulation, CFC-free, to protect product against high ambient conditions
- All-steel cabinet with high-impact epoxy finish for easy cleaning
- Heavy-duty counterbalanced lid with key lock

### Warranty Information

- Domestic and Canada - 13 months parts and labor
- International - 16 months parts

### Accessories

- Optional built-in 6" temperature recorder.
- Optional inventory components, including racks and boxes, may be ordered separately or as a full-load inventory systems

Description	Catalog No.
Inventory Racks ordered separately; Freezer holds four (4) racks	
Rack (holds four 2" Boxes).....	6232
Rack (holds three 3" Boxes).....	6233
2" Fiberboard Boxes (1 doz.).....	5954
3" Fiberboard Boxes (1 doz.).....	5956
6" Temperature Recorder (factory installed)....	6183-6A185
(for "V" and "W" voltages).....	6183-6D185



Temperature Recorder



Inventory Racks

Model No.	Voltage	Amps/Breaker	Cu.Ft./Liters	Interior Dimensions inches/cm	Exterior Dimensions inches/cm	Ship Weight
ULT185-5-A	115V 60Hz1Ø			9.63"H x 12"F-B x 15.25"W	20"H x 24"F-B x 30.75"W	
ULT185-5-V	220V 50Hz1Ø	12.0/15	1.0/28	24.5cm x 30.5cm x 38.7cm	50.8 cm x 61cm x 78.1cm	190 lbs
ULT185-5-W	240V 50Hz1Ø					



**REVCO**  
the future, inside

Revco Technologies, Inc. • 275 Aiken Rd, Asheville NC 28804 USA  
Toll-Free 800-252-7100 • Canada 800-447-3826 • Fax 828-645-3368  
Tel 828-658-2711 • www.revco-sci.com • sales@revco-sci.com  
International Dealers Worldwide

Specifications subject to change without notice.  
©2001 Revco Technologies, Inc. 185R0200 rev0101 Printed in USA.