

## Precision balance Talent TE612



### Standard features.

- High-grade stainless steel weighing pan
- Large, high-contrast LCD
- High level of operating convenience thanks to keys with positive click action and 2 tare keys
- Built-in application programs: counting, weighing in %, net-total formulation; averaging; toggling between 2 units
- Mechanical overload protection
- Bidirectional RS-232C data interface port
- ISO/GLP-compliant logging capability for calibration, adjustment and weighing data in conjunction with the optional YDP03-OCE data printer
- Four digital filter levels for adaptation of the balance to the conditions at the place of use
- Battery-operable (eight 1.5-V Mignon AA batteries) or runs on line current (AC adapter included as part of the equipment supplied)

Readability	0.01 g
Weighing capacity	610 g
Pan size	Ø 116 mm
Response time (avg.)	2 s
Housing (WxDxH)	200x270x70 mm
Reproducibility	<=± 0.01 g
Linearity	<=± 0.02 g



Data printer



SartoConnect

Key accessories	Order no.
Data printer with date, time and statistics functions	YDP03-OCE
Data transfer software SartoConnect with RS-232C standard cable (1.5 m);	YSC01L
Test and calibration weight: 500 g, accuracy class F2	YCW5548-00
External rechargeable battery pack	YRB08Z
Carrying case	YDB01TE
Dust cover	6960TE02