

SPONGEX & LABPAN



Ice buckets - Unbreakable & Light

Unbreakable & Light

Ice buckets and racks

- *made of foamed PVC*
- *no styrofoam*
- *in various colours and types available*

For cooling of experimental products at your working place, we offer the ice buckets SPONGEX and LABPAN as well as suitable racks.

Made of foamed PVC (no styrofoam).

Unbreakable, light and inert against

- *liquid nitrogen*
- *dry ice*
- *salt*
- *acid, base, alcohol, acetone, phenol etc.*

The racks swim on top of liquids and are therefore ideal for use as water bath inserts etc.

They feature a high temperature resistance from - 196 °C to + 93 °C.



SPONGEX

SPONGEX
Ice bucket, round



Including lid.

The racks have a diameter of 220 mm and are available in 2 different sizes:

- for 44 x 1.5/2.0 ml reaction vessels
- for 30 x 14/15 ml FALCON-tubes

LABPAN

LABPAN
Ice bucket, square



Especially developed for loaded racks and vessels of any kind.

LABPAN is available in the following sizes:

MINI (1 litre),
MIDI (4 litres) or
MAXI (9 litres).

The racks are available in 2 different sizes and have the dimensions 235 x 110 mm (w x d):

- for 36 x 1.5/2.0 ml reaction vessels
- for 24 x 14/15 ml FALCON-tubes



Technical Data/Ordering Information

	Capacity litres (below lid)	Total volume litres	Outer Ø at the bottom mm	Outer Ø at the top mm	Outer height mm	Colour	Cat.-No.
SPONGEX	4	5	180	250	180	blue	3.680 052
	4	5	180	250	180	red	3.680 152
	4	5	180	250	180	yellow	3.680 252
Racks	for 44 reaction vessels 1.5/2.0 ml				10	blue	3.680 022
	for 44 reaction vessels 1.5/2.0 ml				10	red	3.680 122
	for 44 reaction vessels 1.5/2.0 ml				10	yellow	3.680 222
	for 30 FALCON-tubes 14/15 ml				17	blue	3.680 042
	for 30 FALCON-tubes 14/15 ml				17	red	3.680 142
	for 30 FALCON-tubes 14/15 ml				17	yellow	3.680 242

	Outer dimensions (w x h x d) mm	Inner dimensions (w x h x d) mm	Height mm	Colour	Cat.-No.	
MINI LABPAN	150 x 80 x 150	100 x 65 x 100	95	blue	3.680 062	
	Capacity 1 litre	150 x 80 x 150	100 x 65 x 100	95	red	3.680 162
		150 x 80 x 150	100 x 65 x 100	95	yellow	3.680 262
MIDI LABPAN	330 x 130 x 230	280 x 115 x 193	125	blue	3.680 072	
	Capacity 4 litres	330 x 130 x 230	280 x 115 x 193	125	red	3.680 172
		330 x 130 x 230	280 x 115 x 193	125	yellow	3.680 272
MAXI LABPAN	410 x 130 x 330	360 x 110 x 250	115	blue	3.680 092	
	Capacity 9 litres	410 x 130 x 330	360 x 110 x 250	115	red	3.680 192
		410 x 130 x 330	360 x 110 x 250	115	yellow	3.680 292
Racks	for 36 reaction vessels 1.5/2.0 ml			10	blue	3.680 302
	for 36 reaction vessels 1.5/2.0 ml			10	red	3.680 402
	for 36 reaction vessels 1.5/2.0 ml			10	yellow	3.680 502
	for 24 FALCON-tubes 14/15 ml			17	blue	3.680 312
	for 24 FALCON-tubes 14/15 ml			17	red	3.680 412
	for 24 FALCON-tubes 14/15 ml			17	yellow	3.680 512

schuett  **biotec.de**

Rudolf-Wissell-Straße 13, D-37079 Göttingen, Germany
 Fon +49 (0) 551/5 04 10-0, Fax +49 (0) 551/5 04 10-99
 info@schuett-biotec.de
 www.schuett-biotec.de