



Laboratory Equipment

Biohazard Shields



Scie-Plas biohazard shields in 6mm acrylic are designed to give laboratory personnel splash protection from biohazardous materials, while providing an unobstructed view of the work surface. Each shield is angled to avoid annoying reflections while an extended base provides stability and a large work area.

Four different models are available: large, small, slim-line and cell counter. The slim-line model is a narrower version of the large model and allows easy access to the work area with both hands. The cell counter is designed for use with small manual and semi-automated haematology counters. Its base slides under the counter for added stability.

ORDERING INFORMATION

Biohazard shield small, 30 x 20 x 37cm	GLE-BS1
Biohazard shield large, 45 x 31 x 52cm	GLE-BS2
Biohazard shield slim-line, 35 x 32 x 52cm	GLE-BS3
Biohazard shield cell counter, 20 x 17 x 26cm	GLE-BS4

Biohazard Flexi-Shields



Scie-Plas biohazard flexi-shields are a lightweight shield mounted on to flexible support arms. These arms are intended to be fixed to the laboratory bench for stability. Once assembled the flexi-shield can be twisted into any shape to provide optimal protection to a technician from splashes when either sitting or standing.

ORDERING INFORMATION

Flexi-shield, 457 x 305mm	GLE-FS1
Flexi-shield, 508 x 406mm	GLE-FS2