

- For the disposal of used pipette tips - additional hole/cover reduces exposure.
- Small and large models
- Waste bags available for easy removal of contents
- Multi-tip boxes feature a slot in the lid, to allow disposal of tips from a multi-channel pipette

TECHNICAL SPECIFICATION

Material	Beta 10mm acrylic	Gamma 12mm lead acrylic (H x W x D)	Dimensions
2L	RPP-B15SC	Ext 150 x 150 x 150mm Int 130 x 130 x 130mm	
2L	RPP-GB15SC	Ext 150 x 150 x 150mm Int 126 x 126 x 126mm	
5.5L	RPP-B32SC	Ext 350 x 150 x 150mm Int 330 x 130 x 130mm	

- Hand protection when pipetting
- For use with Pipetman®
- Convenient to use

Boxes and Pipette Shields

Beta and Gamma tip boxes



Available for either Beta or Gamma radiation protection, these boxes have a hinged lid with a central hole/slot also covered by a smaller hinged lid. The central hole/slot reduces the exposure to radiation during pipetting procedures. Models are available for single or multi-tip disposal. Convenient disposable heavy-duty plastic bags are available.

ORDERING INFORMATION

Tip box-15, Beta (2L capacity) use waste bags RPP-BAG17	RPP-B15SC
Tip box-15, Gamma (2L capacity) use waste bags RPP-BAG17	RPP-GB15SC
Multi-tip box-15, Beta (2L capacity) use waste bags RPP-BAG17	RPP-B15SCMC
Multi-tip box-15, Gamma (2L capacity) use waste bags RPP-BAG17	RPP-GB15SCMC
Tip box-32, Beta (5.5L capacity) use waste bags RPP-BAG35	RPP-B32SC

Pipette Shields



These specifically designed pipette shields fit snugly on the pipette shaft. This protects the hand during the dispensing sequence and enables the user to dispense safely into small and narrow neck containers. Two types of shields are manufactured which fit either Gilson or Biohit pipettes.

ORDERING INFORMATION

Pipette shield, Beta - P20/100, fits Gilson pipettes P20 and P100	RPP-PSG20/100
Pipette shield, Gamma - P20/100, fits Gilson pipettes P20 and P100	RPP-GPSG20/100
Pipette shield, Beta - P200, fits Gilson pipettes P200	RPP-PSG200
Pipette shield, Gamma - P200, fits Gilson pipettes P200	RPP-GPSG200
Pipette shield, Beta - P1000, fits Gilson pipettes P1000	RPP-PSG1000
Pipette shield, Gamma - P1000, fits Gilson pipettes P1000	RPP-GPSG1000
Pipette shield, Beta - PS200, fits Biohit pipettes up to 200µl	RPP-PS200
Pipette shield, Gamma - PS200, fits Biohit pipettes up to 200µl	RPP-GPS200
Pipette shield, Beta - PS1000, fits Biohit pipettes up to 1000 µl	RPP-PS1000
Pipette shield, Gamma - PS1000, fits Biohit pipettes up to 1000µl	RPP-GPS1000
Pipette shield, Beta - PS5000, fits Biohit pipettes up to 1-5ml	RPP-PS5000
Pipette shield, Gamma - PS5000, fits Biohit pipettes up to 1-5ml	RPP-GPS5000