



Desiccator Cabinets

Scie-Plas desiccator cabinets are designed for the storage of anhydrous, hygroscopic material in dust-free and moisture-proof conditions. Each cabinet is provided with vented shelves that maximise air flow and drying efficiency. Although available in 3 standard sizes, the dimensions of these cabinets can be customised to meet your own personal requirements - just contact sales@scie-plas.co.uk for a quotation.

BENEFITS INCLUDE

- **Available in 3 sizes** - large (DESC-LG), medium (DESC-MD) and small (DESC-SM)
- **Robust construction** - from transparent 10mm acrylic
- **Gas tight silicone seal** - recessed into each cabinet door to maintain air- and moisture-free conditions
- **Vented shelves (DESC-SH)** - maximise air flow and drying efficiency
- **Non-vented, lipped desiccant tray (DESC-DS)** - allows desiccant fluid or powder to be stored in the bottom of the cabinet

ORDERING INFORMATION

Power Supply

	Part No.
Small desiccator with 2 vented shelves and 1 desiccant tray	DESC-SM
Medium desiccator with 3 vented shelves and 1 desiccant tray	DESC-MD
Large desiccator with 4 vented shelves and 1 desiccant tray	DESC-LG

Replacement Parts & Accessories

1 x vented shelf	DESC-SH
1 x non-vented, lipped desiccant tray	DESC-D5

TECHNICAL SPECIFICATION

	DESC-SM	DESC-MD	DESC-LG	DESC-SH	DESC-DS
Dimensions (W x H x D)	305 x 178 x 305mm	305 x 305 x 305mm	305 x 457 x 305mm	290 x 17 x 285mm	280 x 25 x 280mm
Internal Dimensions (W x H x D)	285 x 158 x 285mm	285 x 285 x 285mm	285 x 437 x 285mm	-	-
Number of Shelves, Including Desiccant Tray	3	4	5	-	-
Number of Shelf Slots	4	8	12	-	-
Weight	6 kg	7.8 kg	8.6 kg	0.2 kg	0.3 kg



Glass Plate Racks

These sturdy racks are designed for safe drying and storage of glass plates. The small rack can hold up to twenty 2mm thick plates, while the larger rack can take up to ten 4mm thick glass plates.

ORDERING INFORMATION

Mini glass plate rack for twenty 2mm plates	GLE-MGR
Large glass plate rack for ten 4mm plates	GLE-LGR