

Microtube Racks, Combi-Racks, Mini-Racks

Combi-Racks

Each combi-rack is configured in 6 rows of 12 holes. Row 1 is for 1.5ml microtubes and rows 2 - 6 are for 0.5ml tubes. This layout is ideal for those using PCR® thermal cyclers which have heating blocks usually accommodating 60 tubes.

ORDERING INFORMATION

Combi-rack - yellow	GLE-TR6Y
Combi-rack - red	GLE-TR6R
Combi-rack - green	GLE-TR6G
Combi-rack - blue	GLE-TR6B
Combi-rack - pk/4 one of each colour	GLE-TR6P

Microtube Racks

Scie-Plas general purpose racks are ideal for everyday usage of 0.5ml and 1.5ml microcentrifuge tubes as well as screw cap tubes and cryovials. They are made from durable PVC and are suitable for autoclaving. Available in four different colours for easy user identification on the bench or in a refrigerator.

ORDERING INFORMATION

Rack for 100 x 1.5ml tube - yellow	GLE-TR3Y
Rack for 100 x 1.5ml tube - red	GLE-TR3R
Rack for 100 x 1.5ml tube - green	GLE-TR3G
Rack for 100 x 1.5ml tube - blue	GLE-TR3B
Rack for 100 x 1.5ml tube - pk/4 one of each colour	GLE-TR3P
Rack for 50 x 1.5ml tube - yellow	GLE-TR2Y
Rack for 50 x 1.5ml tube - red	GLE-TR2R
Rack for 50 x 1.5ml tube - green	GLE-TR2G
Rack for 50 x 1.5ml tube - blue	GLE-TR2B
Rack for 50 x 1.5ml tube - pk/4 one of each colour	GLE-TR2P
Rack for 60 x 0.5ml tube - yellow	GLE-TR1Y
Rack for 60 x 0.5ml tube - red	GLE-TR1R
Rack for 60 x 0.5ml tube - green	GLE-TR1G
Rack for 60 x 0.5ml tube - blue	GLE-TR1B
Rack for 60 x 0.5ml tube - pk/4 one of each colour	GLE-TR1P

Mini-Racks

Scie-Plas mini-racks are small convenient racks for those working with just a few 0.5ml or 1.5ml microtubes. Having no base, they are ideal for use with ice buckets, allowing direct contact with ice.

ORDERING INFORMATION

Clear mini-rack 12 x 1.5ml tubes	GLE-TR4
Clear mini-rack 16 x 0.5ml tubes	GLE-TR5

Microtube Rack Storage

These acrylic stores for popular 80-well microtube racks are designed for the convenient short term storage of large numbers of microtubes. They provide a convenient solution to sample organisation on the workbench top or in a refrigerator. Each compartment is generously sized to enable easy removal or replacement of the racks. Both vertical and horizontal options accommodate up to eight racks, allowing up to 640 x 1.5/2ml microtubes or similar tubes to be stored locally.

ORDERING INFORMATION

8-place vertical storage rack for 80-well microtube racks	GLE-VSR
8-place horizontal storage rack for 80-well microtube racks	GLE-HSR

