

Workstation / Storage Bin

This acrylic unit is an ideal method of storing and dispensing laboratory disposables such as pipette tips, microtubes, medical wipes, disposable gloves and sealing films, eg. Parafilm®, in one central location. The top compartments feature hinged lids to keep tips and tubes free from dust and contamination. The unit can either sit on the bench or be mounted on a wall.

ORDERING INFORMATION

Workstation storage bin

GLE-WSB



Medical Wipes Dispensers

By locating these dispensers on a bench or wall surface, laboratory wipes can be single-handedly removed from their original packaging. With internal dimensions of 3.5 x 13 x 11cm (small) and 5.5 x 32 x 16.5cm (large), most manufacturers' products can be accommodated.

ORDERING INFORMATION

Medical wipes box holder - small

GLE-KBS

Medical wipes holder - large

GLE-KBL



Safety Spectacles Dispenser

The dispenser is fabricated from 3mm acrylic and keeps safety spectacles dust-free, visible and available to laboratory personnel and visitors. Spectacles are removed from the bottom opening and replaced through the top afterwards.

ORDERING INFORMATION

Safety spectacle dispenser

GLE-GB



Petri Dish Dispenser

The Scie-Plas petri dish dispenser can sit on a bench or be fixed to a wall surface. They allow easy one-handed access to petri dishes up to 98mm diameter. The two compartments will hold up to 24 standard 90mm dishes

ORDERING INFORMATION

Petri dish dispenser

GLE-PDD



Serological Pipette Rack

These racks are for bench top storage of manual glass or plastic serological pipettes of 0.1ml to 25ml capacity. They have four angled compartments measuring 8 x 8 x 30cm. being fabricated from 5mm clear acrylic, stock levels of pre-packed pipettes can be constantly monitored. Additionally, each compartment can accept standard pipette cans and allow one-handed retrieval of pipettes when working in a biological hood.

ORDERING INFORMATION

Pipette rack - four shelves

GLE-PR



Parafilm® Sealing Film Dispenser

The Scie-Plas sealing film dispenser is of a robust clear acrylic construction to allow easy monitoring of remaining film. A sharp blade is safely located under the hinged lid for user protection. The dispenser can accommodate 1 x 100mm or 2 x 50mm wide rolls of Parafilm® or similar sealing films.

ORDERING INFORMATION

Parafilm® dispenser 1 x 100mm or 2 x 50mm

GLE-PD



Pipette Stands and Holders

Pipette Stands



Designed to keep pipettes off potentially contaminated surfaces and reduce bench clutter, Scie-Plas pipette stands are available in 3 or 4-position formats. In attractive coloured acrylic, they have been designed to accept most brands of pipettes (slot width 25mm).

The multi-channel stand allows any pipette (multi or single) with a hooked finger-rest to be simply placed on the crossbar between the two upright base plates. This stand's capacity is limited only by width and will accommodate up to 3 x standard 8-channel pipettes.



ORDERING INFORMATION

3-pipette stand	GLE-PS3
4-pipette stand	GLE-PS4
Multi-channel pipette stand	GLE-PSMC



Pipette Filler Holder

The Scie-Plas pipette filler holder is designed to hold motorised pipetting aids. The 25mm locating hole will accept the pipette holder shaft from most models of 8mm acrylic and with a solid base of 14 x 14cm, the holder is extremely stable and durable.

ORDERING INFORMATION

Pipette filler holder	GLE-PFH
-----------------------	---------



Draught Protection Cabinets

Scie-Plas draught protection cabinets are simple chambers designed to protect sensitive instruments such as microscopes or balances from the local environment. The balance cabinets feature hinged lids for easy access to the weighing pan. Fabricated from ultra-light 3mm transparent acrylic, the read out is easily visible when the lid is closed.

ORDERING INFORMATION

Balance/microscope protector 250 x 450 x 480cm	GLE-MC1
Balance cover, hinged lid and hole 220 x 330 x 160cm	GLE-BC1
Balance cover, hinged lid and hole 260 x 290 x 160cm	GLE-BC2
Balance cover, lift-off lid and 2 handles 330 x 300 x 310cm	GLE-BC3
Balance cover, lift-off lid and 2 handles 410 x 410 x 550cm	GLE-BC4



Coloured Tape Dispensers

Scie-Plas tape dispensers are ideal for easy one-handed access to popular coloured identification tapes. A serration on the dispensing edge allows the tape to be easily torn off. Each dispenser will take up to 12 rolls of the most popular 12.7mm wide tapes. Of smoked acrylic construction, two models are available for tapes of differing core diameters.

ORDERING INFORMATION

Small tape dispenser, tape core diameter 26mm	DISP26MM
Large tape dispenser, tape core diameter 76mm	DISP76MM