

# Wireless Alarm System



## Data Loggers (8 bit resolution) with Digital Display

Logger Type (A) **Range -200C to +100C**

Suitable for Liquid Nitrogen Vessels and U.L.T Freezers with digital display (8 bit logger, temperature steps 1.17 over range with accuracy +/- 0.3C)

Logger Type (B) **Range -30C to +50C**

Suitable for all Refrigeration ,Incubators, Bloodbanks etc.

(8 bit logger, temperature steps 0.31 over range with accuracy +/- 0.2C)

**Replaces Chart Recorders**

supplied with full DIY instructions



## Data Loggers (12 bit resolution)

Logger Type (C) **Range -40C to +125C**

Suitable for high resolution reading requirements (all refrigerated lab.equipment) (12 bit logger, temperature steps 0.04 over range with accuracy +/- 0.2oc)

**On site Demonstration on request**

humidity, pressure and digital signal loggers also available

**CRYOTECHNICS**

*the science of refrigeration.....*