

TENAK EL51LT FREEZER



Technical specifications

Model	EL51LT
Article no	TE21128
Gross volume, lt.	416
Gross volume, cu.ft.	14,9
Insulation thickness	100
No. of racks	30 (27 + 3 on bench)
Temperature ambient 10°C to 30°C	-10°C to -45°C
Time for reaching -45°C at start up, empty	3 hours
Temperature stability on filled freezer	± 0,5°C
Energy consumption	0,48 kWh
Noise level	<48 dB
Power supply	220/230V - 50 Hz
Electric circuit Fuse	10A
Power input	3 – 4 A
Shipping dimensions (HxDxW)	980 x 770 x 1760 mm
Shipping volumen	1,26 cbm
Shipping weight	81 kg
Inner dimensions (HxDxW)	650 mm x 450 mm x 1500/1300 mm
Maximum rack height	640 mm
Outer dimensions without castors (HxDxW)	850mm x 655 mm x 1705 mm
Outer height with castors	865 mm
Net weight	77 kg
Tariff number	84 18 30 90

• Article no. TE21128


The EL51LT stores up to 31500 2 ml cryovials by using:

TE21450	3 ea	Racks for 6 boxes 50 mm
TE21454	27 ea	Racks for 11 boxes 50 mm
TE22005	315 ea	White cardboard box 50 mm
TE22120	315 ea	10x10 divider for 50 mm box
TE78020	315 ea	2 ml cryovials 100/pack
TE78106	63 ea	Ass. colour inserts for vials 500/pack

See "COOL AND SAFE STORAGE" for prices and more options.

Inside measures

A =	650 mm
B =	450 mm
C =	1500 mm
D =	400 mm
E =	200 mm
F =	250 mm
G =	1300 mm



Features

- Acoustic Alarm
- Large temperature display
- Easy operated Dixell controller for temperature setting between -10°C and -45°C
- Electronic thermostat
- Heavy duty castors for easy moving (optional, not included)
- Handle with key lock
- Defrost water drain hole

Additional information

- CFC free refrigerant
- One Stage Danfoss compressor (Custom design) ensures high level of reliability
- Space saving design
- Double seal lid gaskets
- Is manufactured according to ISO9000 and is CE marked
- A high standard test procedure on each product ensures quality freezers
- Supplied with Schuko plug as standard

TENAK