

WOLF LABORATORIES

Tel: 01759 301142 Fax: 01759 301143 Email: sales@wolflabs.co.uk
Website: www.wolflabs.co.uk

Mineralight UV Display Lamps



Mineralight® display lamps use two 25-watt high performance tubes providing hours of UV illumination. Each lamp comes with adjustable handles and mounting bracket for hanging

- Two 25-watt UV tubes for bright illumination
- Available with 254nm, 302nm or 365nm UV
- Use the brackets included with the lamp to attach the lamp to the ceiling or other structure for overhead UV illumination

Specifications

UV Source: 254nm, 302nm or 365nm UV

Watts: Twenty-five

Dimensions: 18.5L x 4.75W x 5H in. (470 x 114 x 89mm)

Weight: 10 lb. (4.6 kg)

Applications

- Fluorescence
- Sterilization
- Inspection