



Ind3530 Ultrasonic Bath Complete Spec Sheet

SPECIFICATIONS

DETAILS

CONSTRUCTION CASE	RIGID STAINLESS STEEL
CONSTRUCTION TANK	STAINLESS STEEL WELDED 2MM
OPERATING VOLTAGE	220/240 VAC 1PH 50HZ
POWER COMSUMPTION ULTRA.	400 VA RMS
POWER COMSUMPTION HEATER	800 VA (2 X 400 HEATERS)
DRY RUN PROTECTION	STAINLESS STEEL FLOAT SWITCH
TRANSDUCERS METAL BONDED	8 X 50MM PILLARS
POWER OUTPUT INTO FLUID	16 WATT / LITRE
THERMOSTAT CONTROL	30 - 75 oC
TIMER RANGE (ON GENERATOR)	0 - 15 MINS
TANK DIMENSIONS	350 X 300 X 300 DEEP
MAX CAPACITY	31.50 LITRES
WORKING CAPACITY	25.00 LITRES
DRAIN OUTLET	3/4" BSP VALVE
OPERATING FREQUENCY	30 - 40 KHZ
EFFICIENCY	83%
OVERALL DIMENSIONS	500 X 450 X 400 HIGH
NET WEIGHT	35.0 KGS
GROSS WEIGHT FOR EXPORT WITH LID, BASKET AND PACKING	55.0 KGS

ACCESSORIES AVAILABLE:

STAINLESS LID, BASKET, CLEANING FLUID, + STAINLESS STEEL STAND