



# Ultrawave Ltd

## Ind500 Ultrasonic Bath Complete Spec Sheet

### SPECIFICATIONS

### DETAILS



CONSTRUCTION CASE	RIGID STAINLESS STEEL
CONSTRUCTION TANK	STAINLESS STEEL WELDED 2MM
OPERATING VOLTAGE	220/240 VAC 1PH 50HZ
POWER COMSUMPTION ULTRA.	150 VA RMS
POWER COMSUMPTION HEATER	150 VA
HEATER RANGE DIGITAL RANGE	0 - 80 oC
TRANSDUCERS METAL BONDED	3 X 50MM PILLARS
POWER OUTPUT P-P WATTS	300 WATTS
TIMER RANGE DIGITAL DISPLAY	0 - 59 MINS
TANK DIMENSIONS	300 X 150 X 150 DEEP
MAX CAPACITY	5.00 LITRES
WORKING CAPACITY	4.40 LITRES
OPERATING FREQUENCY	30 - 40 KHZ
MODULATION	50 - 100 HZ
OVERALL DIMENSIONS	330 X 180 X 285 HIGH
NET WEIGHT	6.6 KGS
GROSS WEIGHT FOR EXPORT WITH LID, BASKET AND PACKING	9.7 KGS
DRAIN OUTLET	1/2" BSP VALVE

**ACCESSORIES AVAILABLE:**

STAINLESS LID, BASKET, CLEANING FLUID

Tel: 0845 330 4236 / Fax: 0845 330 4231  
W: [www.ultrawave.co.uk](http://www.ultrawave.co.uk) / E: [sales@ultrawave.co.uk](mailto:sales@ultrawave.co.uk)