

Sonicating Water Baths



Unit Type	U50
Litre Capacity	0.5
Transducers	1 Disc type
P-P Watts	100
Heat	None
Nett Weight	1.4kg
Internal mm	150x85x65
Drain	None
Timer	0-15 mins Analogue



Unit Type	U300H
Litre Capacity	2.75
Transducers	2 Disc type
P-P Watts	200
Heat	30-70°C
Nett Weight	4.5 Kg
Internal mm	240x140x100
Drain	None
Timer	0-15 mins Analogue



Unit Type	U500D
Litre Capacity	4.75
Transducers	2 Pillar type
P-P Watts	280
Heat	0-80°C Digital
Nett Weight	5.6 Kg
Internal mm	300x150x150
Drain	Yes
Timer	0-59 mins Digital



Unit Type	U2800D
Litre Capacity	24.0
Transducers	8 Pillar type
P-P Watts	800
Heat	0-80°C Digital
Nett Weight	18.5 Kg
Internal mm	500x300x200
Drain	Yes
Timer	0-59 mins Digital