

replacement ultraviolet tubes

UVItec can supply a full range of ultraviolet tubes to fit most available models of UV transilluminator and UV lamp, plus other UV devices. We specialise in rapid delivery from extensive stocks. We also offer discounts for purchase of multiple (5+) tubes.

In addition to those listed in the table, we can also supply larger UV tubes of up to 40W including short wave (254nm) germicidal tubes.

Please contact us for details, or see our web site for a complete list.

Ordering

Model No.	Output (Watts)	Wavelength
T4L	4W	365
T4C	4W	254
T4LN	4W	365 BLB
T6L	6W	365
T6C	6W	254
T6M	6W	312
T6LN	6W	365 BLB
T6LC	6W	365/254
T8L	8W	365
T8C	8W	254
T8M	8W	312
T8LN	8W	365 BLB
T15L	15W	365
T15C	15W	254
T15M	15W	312
T15LN	15W	365 BLB



ultraviolet safetywear

Safetywear to protect the eyes and face from harmful ultraviolet rays is recommended to all operators of ultraviolet equipment and patients undergoing PUVA treatment. Cover spectacles and Face Shields are specially produced for UVItec by Bollé, the leading European eye safetywear manufacturer.

- Comfortable to wear
- Lightweight, unobtrusive design
- Guaranteed to completely block ultraviolet radiation below 400nm
- CE marked

Ordering


SS350	UV safety cover-specs
SG550	UV safety goggles
SF850	UV face shield

uvisafe



uvigene

High quality thermal cycler for 0.2 or 0.5ml microtubes



uvipage

uvisub

Professional vertical and horizontal electrophoresis equipment

for more information on the above products

visit our website at www.uvitec.co.uk

call us **+44 (0)1223 568060** email uvi@uvitec.co.uk