

Pipettes - micro / bulbs

Pipettes, micro-sampling

03



IntraMARK micropipettes, BLAUBRAND®, with ring mark are according to DIN ISO 7550, and are conformity certified for volumes $\geq 5\mu\text{L}$: Accuracy $\geq \pm 0.25\%$ Reproducibility $\geq 0.5\%$ (size $5\mu\text{L}$: Accuracy $\geq \pm 0.3\%$; Reproducibility $\geq 0.6\%$) Fire polished surfaces inside and outside are particularly smooth, the aspiration end is also fire polished to minimise pipette adapter wear. Calibrated to contain (TC, In) filled from one end to a ring mark. Colour code to ISO for clear identification. Length 127mm (size $200\mu\text{L}$: 140mm). Pack quantity 250. For use with PMK-850-H below.

Catalogue No **Alt. No** **Capacity, μL** **Colour code**

PMK-800-010Q	708707	1/2/3/4/5 [†]	White
PMK-800-030K	708709	10	Orange
PMK-800-050E	708718	20	Black
PMK-800-070V	708722	25	2 x White
PMK-800-090P	708733	50	Green
PMK-800-110M	708744	100	Blue
PMK-800-130G	708757	200	Red

[†] Conformity certified for $5\mu\text{L}$ mark.

Micro-pipettor aid

03



Indispensable accessory for sampling with disposable micropipettes with ring mark and many small volume pipettes up to 1mL (e.g. blood diluting and blood sugar pipettes). The micro reduces the hazards of infection and is autoclavable at 121°C . The integral ejection device allows the disposal of contaminated pipettes without touching them. The micro is extremely light and convenient.

Catalogue No
PMK-850-H

Pipettes, bulbs, PVC, pear shape

02



Pack quantity 10.

Catalogue No	Bulb diameter, mm	Neck bore, mm	Capacity, mL
PMR-200-010R	34	4	15
PMR-200-030L	39	4	22
PMR-200-050F	43	4	30
PMR-200-070W	50	4	42

Pipettes, bulbs, rubber teats

02



Catalogue No	Neck bore, mm	Nominal capacity, mL	Pack qty
PMR-210-010K	5	1.5	20
PMR-210-030E	5	2	20
PMR-210-050V	6	3	20
PMR-210-070P	6	5	20

Catalogue No	Description	Pack qty
PMR-220-G	Mixed pack of 5 of each of the above	20

Pipettes, fillers, Pi-pump

02



An almost indestructible, acid resistant pipette controller, easily dismantled for cleaning. Suitable for filling and dispensing from all standard pipettes from micro to 25mL.

Catalogue No	Colour	For pipette sizes, mL
PMR-336-005X	Yellow	0.2
PMR-336-010V	Blue	2
PMR-336-030P	Green	10
PMR-336-050J	Red	25