

Fibreboard Freezer Boxes



Economical cryovial boxes fit all liquid nitrogen and mechanical square or upright storage racks. Select between two box sizes (lids included) as well as 100-, 81-, 64- and 49-cell partitions that fit in either box. All the boxes and the partitions sell in packs of 12.

WO-A1	Fibreboard freezer boxes 130x130x50mm - pack of 12
WO-A2	Fibreboard freezer boxes 130x130x75mm - pack of 12
WO-A6	Fibreboard freezer boxes 130x130x95mm - pack of 12
WO-A1/49DIVS	16mm 49 cell dividers For A1,A2 & A6 - pack of 12
WO-A1/64DIVS	14mm 64 cell dividers For A1,A2 & A6 - pack of 12
WO-A1/81DIVS	13mm 81 cell dividers For A1,A2 & A6 - pack of 12
WO-A1/100DIVS	12mm 100 cell dividers For A1,A2 & A6 - pack of 12