



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax:01759 301143

sales@wolflabs.co.uk



If the literature below doesn't answer all your questions, click these links to:

[Request details on an alternative product of the same type](#)

[Request details on accessories](#)

[Request details on other product categories](#)



High Accuracy Probe Thermometer



High Accuracy Probe Thermometer

ST-35

Features:

- Measures in °C
 - Fast response and accuracy as required by food hygiene regulations
 - Water-resistant instrument with Lumberg connector for improved prevention of food and liquid ingress and more secure probe connections
 - General purpose 100mm probe, 3.3mm diameter with 0.7m coiled cable
 - 1 year warranty on thermometer
 - Large LCD display
 - Auto off after 3 mins
- **Supplied with general purpose probe (ST35-GP), rubber boot, 5 coloured thumb tabs and battery**

• **Optional Calibration Certificate available**

CAL-STSYS Calibration Certificate - Issued at 4 points: 25°C, 50°C, 75°C & 100°C. (Additional points are available if required. These must be specified at the time of order. Please call for price and details.)

Recommended Temperature Probes

ST35-AP Flexible Air Probe. Thermistor sensor with straight cable. Lumberg Connector. Time Constant 1.9 secs. Max. Tip Temp 120°C/248°F

ST35-GP General Purpose Probe. Thermistor sensor with coiled cable 100mm, 3.3mm diameter, 0.7m straight cable. Lumberg Connector. Time constant 1.3 secs. Max. Tip Temp. 120°C/248°F

High Accuracy Probe Thermometer

Technical Specifications

Range	-40° to +110°C
Resolution	0.1°
Accuracy	± 0.5°C (0°C to 70°C), ± 1°C remaining range
Power	1 x 9V PP3 battery
Dimensions	70 x 165 x 35mm
Weight	172g