



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax:01759 301143

sales@wolflabs.co.uk



If the literature below doesn't answer all your questions, click these links to:

[Request details on an alternative product of the same type](#)

[Request details on accessories](#)

[Request details on other product categories](#)



Probe Thermometer



Probe Thermometer

P-300

Features:

- Measures in °C or °F
- Precision NTC thermistor sensor for fast accurate readings
- Dust-protected and splash-proof to IP 54 classification (IEC standard 60529)
- Suitable for food temperature measuring to HACCP regulations
- Data hold retains the displayed reading
- Auto power off after 10 minutes of no activity
- Electro-luminescent backlight for use in low light environments
- Impact-resistant ABS case with flip-out desk stand to use free standing
- Splash-proof membrane keypad
- 120mm stainless steel probe connected via 1 metre cable
- 22mm LCD display with function indicators

• Supplied with probe & 2 batteries

• Optional Calibration Certificate available

CAL-ST Calibration Certificate - Issued at 4 points: 25°C, 50°C, 75°C & 100°C. (Additional points are available if required. These must be specified at the time of order. Please call for price and details.)

Probe Thermometer

Technical Specifications

Range	-40° to 200°C	-4° to 392°F
Resolution	0.1°C	0.1°F
Accuracy	±0.5°C (0° to 100°C),	±2° (above 100°C)
Sampling	every 0.5 second	
Power	2 x 1.5V AA batteries ALR-6	
Dimensions	260 x 96 x 25mm	
Weight	187g	

Dual Sensor Memory Indoor Outdoor Thermometer



Dual Sensor Memory Indoor Outdoor Thermometer

TH-968

Features:

- Measures °C or °F
- Easy to operate
- Min/max memory with separate reset buttons for independent reset for internal and external sensors
- Water-resistant external sensor with 3m connection lead
- Separate windows with 15mm LCD displays for each sensor with IN and OUT indicators
- Impact-resistant ABS case with flip-out desk stand and wall mounting key-ways

· Supplied with battery

· **CAL-TH Optional Calibration Certificate available** Issued at 2 points for external sensor: -15°C & 20°C and 2 points for internal sensor: 0°C & 20°C. Additional points available if required. Please call for details

Dual Sensor Memory Thermometer

Technical Specifications

Range	-50° to +70°C, -58° to +158°F (external sensor) -10° to +50°C, 14° to +122°F (internal sensor)	
Resolution	0.1°C	0.1°F
Accuracy	±2°C	±3.6°F
Power	1 x 1.5V AA battery	
Dimensions	135 x 75 x 18mm	
Weight	117g	