



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax:01759 301143

sales@wolflabs.co.uk




If the literature below doesn't answer all your questions, click these links to:

[Request details on an alternative product of the same type](#)

[Request details on accessories](#)

[Request details on other product categories](#)





Bioreader 5000 Pro-V Range

Latest, Generation 5 Bioreader® software level 10

High Resolution true color Elispot + Plaque + CFU Reader for Infection biology, Genetic Toxicology, 10cm Petri dish, Antibiotic potency

With highest resolution digital Megapixel Fire Wire Camera with SONY Chip

Up to 3 colors may be counted in each well, unlimited number may be selected

Versatile Reader for

Elispot assays

Plaque Assays

CFU

- Easy change between different types of plates.
- Confirms CE/ EMV
- Automated focus, key locked
- Software controlled illumination
- Easy to use, professional and GLP user interface
- Triple Lens changer
- Digital camera, 2000*1500 resolution
- Flat Screen Monitor 24 inches 1920*1600 Pixels
- MS-Office (Word/Excel/Outlook)
- DELL Dual Core 2*320 GB Hard disk RAID

Model [Part No.]	Purpose	User interface	Zoom / lens changer	Camera	DELL PC	Dell Monitor
Pro- V alpha [1471]	Elispot/Plaque 96 + 12/24 wells + Zoom Lens system	'Ease of use'	Double changer incl. 2 lenses	1024 * 768	Dual core 160 GB HDD	20 '' 1600 * 1200
Pro- V beta [1674]	Elispot/CFU/Plaques 96 + 12/24 + 6 wells	'Ease of use' / 'professional'	Triple changer incl. 3 lenses	1280 * 960	DUAL CORE 320 GB HDD	20 '' 1600 * 1200
Pro- V gamma [1733]	Infections biology, Genetic Toxicology, Plaque, 10 cm Petridish and 96 /12/24 well plates. Elispot, Petri dishes, Antibiotic potency	'Ease of use' / 'professional' / 'GLP'	Triple changer incl. 3 lenses	2000 * 1500	Quad core 2* 320 GB HDD	24'' 1920 * 1600