



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax:01759 301143

sales@wolflabs.co.uk



If the literature below doesn't answer all your questions, click these links to:

[Request details on an alternative product of the same type](#)

[Request details on accessories](#)

[Request details on other product categories](#)



Biohit Proline® Biocontrol – decontamination solution

Biohit Proline Biocontrol is a safe, simple solution for decontamination in the laboratory. It is specifically designed for decontaminating laboratory instruments.

The components of the pipette, which most often become contaminated, whether by incorrect pipetting techniques or by accident, are in the tip cone area.

Daily safety maintenance

Lightly spray the outside of the pipette and wipe dry with a clean cloth. Clean surfaces, such as incubators and benches with no noticeable contamination such as blood stains may be dampened using the Proline Biocontrol, left on for 2 minutes and then wiped clean with an absorbent cloth.

Effective against:

• Viruses

Hepatitis A, B, C, D and E, HIV, Parvovirus, Newcastle Disease

• Bacteria

Pseudomonas sp., coli, Staph aureus, Proteus mirabilis, Proteus vulgaris, Streptococcus sp., Bacillus sp., Salmonellae, Klebsiella sp., Enterobacter sp., Serratia sp., Mycobacterium sp., Listeria sp., Legionella pneumophila

• Fungi

A. Niger, Penicillium sp., T. Mentagrophytes, Candida albicans

Cat. No.	Qty/Unit	Pack Format
724004	5 litre	Economic Large Container

Biohit Proline® Biocontrol is:

- Non-toxic and non-irritating to humans and animals
- Ready to use
- Lethal to a broad spectrum of viruses, bacteria and fungi
- BS 6471 Approved
- Non-corrosive to metals
- Non-fuming
- Odourless and non-staining
- Non-mutagenic
- Alcohol, aldehyde, phenol and chlorine free
- Effective at removing odours
- Environmentally friendly (biodegradable)
- Stable for long periods



See complete pipette decontamination procedure on page 83

