



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax:01759 301143

sales@wolflabs.co.uk



If the literature below doesn't answer all your questions, click these links to:

[Request details on an alternative product of the same type](#)

[Request details on accessories](#)

[Request details on other product categories](#)





Ideal for laboratory use, featuring accurate temperature control using a digital temperature control with electronic display. Fitted with castors as standard.







Climate Class N

Features

- Integral Lock
- Castors
- Digital temperature controller
- Temperature failure alarm
- Easy clean internal rounded corners



Technical Specifications

Name	Dimensions W x D x H	Internal Dimensions W x D x H	Capacity CU FT (L)	Temp Range
 GLE10	620x400x925	480x250x380	1.9 (54)	-20 / -45°C
 GLE20	980x680x925	800x440x710	7.9 (215)	-20 / -45°C
 GLE30	1280x680x925	1100x440x710	11.2 (310)	-20 / -45°C
 GLE40	1500x680x925	1320x440x710	13.4 (380)	-20 / -45°C
 GLK20	980x690x925	800x440x710	7.6 (215)	-25 / -60°C
 GLK30	1280x690x925	1100x440x710	11.0 (310)	-25 / -60°C

GLE10	Low Temp Freezer	White
GLE20	Low Temp Freezer	White
GLE30	Low Temp Freezer	White
GLE40	Low Temp Freezer	White
GLK20	Ultra Low Temp Freezer	White
GLK30	Ultra Low Temp Freezer	White

Accessories

Name	Description	Product
Basket GLE	Basket	GLE10
Basket GLE	Basket	GLE20
Basket GLE	Basket	GLE30
Basket GLE	Basket	GLE40