



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax:01759 301143

sales@wolflabs.co.uk



If the literature below doesn't answer all your questions, click these links to:

[Request details on an alternative product of the same type](#)

[Request details on accessories](#)

[Request details on other product categories](#)



Microtube Freezer Racks



Inexpensive freezer storage

Lightweight 1.5 to 2.0 ml microtube racks are ideal for long-term studies or tube storage. Two-piece racks are constructed of expanded polystyrene (EPS) foam. Extra spacing between wells allows easy gripping with fingertips. Racks stack securely with nesting features on lid and base.
Pack of 10.

Cat No.	Wells	L x W x H, mm	Units
ZXZHS4283	50	210 x 110 x 71	Pack of 10
ZXZHS4282	100	336 x 95 x 73	Pack of 10