



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax:01759 301143

sales@wolflabs.co.uk



If the literature below doesn't answer all your questions, click these links to:

[Request details on an alternative product of the same type](#)

[Request details on accessories](#)

[Request details on other product categories](#)



WOLF LABORATORIES

Tel: 01759 301142 Fax: 01759 301143 Email: sales@wolflabs.co.uk
Website: www.wolflabs.co.uk

Microtube Storage Rack 100 Well



Durable storage for microtubes and cryogenic vials

Rugged polypropylene box has durable three-point hinges and clasp for secure closure, moulded grid lines on the lid, imprinted coordinates on the bottom of base, and imprinted and raised coordinates on each well

Holds 1.5 to 2.0 ml tubes and vials

Autoclavable

Measures 141 x 151 x 57 mm

Cat No.	Colour	Units
HS1064	Natural	Pack of 5
HS1065	Blue, green, purple, yellow and orange	Pack of 5