



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax:01759 301143

sales@wolflabs.co.uk



If the literature below doesn't answer all your questions, click these links to:

[Request details on an alternative product of the same type](#)

[Request details on accessories](#)

[Request details on other product categories](#)



WOLF LABORATORIES

Tel: 01759 301142 Fax: 01759 301143 Email: sales@wolflabs.co.uk
Website: www.wolflabs.co.uk

PCR Rack 96 Well



Removable hinged lid lies flat. Perfect for pre- or post-PCR sampling. Sturdy polypropylene rack accommodates individual 0.2 ml tubes or strips of 8 or 12 tubes. Tube wells are easily identified with embossed numbers and letters. Rack measures 130 x 98 x 33 mm. Autoclavable.

Cat No.	Description	Units
HS23461A	96-Well PCR Rack, blue, green, pink, yellow and orange	Pack of 5
HS23461B	96-Well PCR Rack, Natural	Pack of 5