



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax:01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**If the literature below doesn't answer all your questions, click these links to:**

[Request details on an alternative product of the same type](#)

[Request details on accessories](#)

[Request details on other product categories](#)



# WOLF LABORATORIES

Tel: 01759 301142 Fax: 01759 301143 Email: sales@wolflabs.co.uk  
Website: www.wolflabs.co.uk

## Rota-Filler® 3000



### Patented, adjustable nozzle rotates to three positions for maximum comfort and convenience

The Rota-Filler 3000 combines a powerful motor with a lightweight, ergonomic design; a reflection of our commitment to the comfort of the end user. Features such as textured grip, hand-neutral design for right- and left-handed users, concave triggers, and convenient controls ensure the filler provides many hours of stress-free filling and dispensing.

Filling and dispensing speed is easily set with a single switch. Force liquid out of your pipette by setting the unit to low- or high-speed blow-out mode, or flip to gravity mode for use with “to contain” pipettes and dispense liquid without the motor. Filler accommodates glass and plastic pipettes ranging in volume from 1 to 100 mL. The filter in the polypropylene nozzle prevents contamination.

Pipetting speed is easily set. Slide the speed switch between two operating modes—blow-out or gravity. Select the low or high range in blow-out mode to dispense the remaining liquid in the pipette. Flip to gravity mode for use with “to contain” pipettes and dispense the liquid without the motor.

A red low-battery indicator signals when the batteries need to be recharged. Filler can be recharged while in use, or while resting in either the dual-position desktop or magnetic wall-mountable holder (both included). Includes two replacement 0.45 µm filters.

### Specifications

**Power:** two NiMH batteries, replaceable

**Filter:** hydrophobic 0.45 µm

**Pipettes:** 1 to 100 mL, glass and plastic

**Pipette nozzle:** polypropylene, autoclavable

**Dimensions:** 170 x 47 x 143 mm

Cat No.	Description	Plug Type	Units
HS3000RFBK	Rota-Filler 3000, 220/230 VAC	U.K.	Each
HS5000RFFIL	Replacement Filters, Pack of 5	N/A	Pack of 5