



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax:01759 301143

sales@wolflabs.co.uk



If the literature below doesn't answer all your questions, click these links to:

[Request details on an alternative product of the same type](#)

[Request details on accessories](#)

[Request details on other product categories](#)



PHARMACY SOLID DOOR. Our pharmacy refrigerators have been designed specifically for pharmacy use and offer a safe and reliable environment for the storage of vaccines and drugs. They adhere to the guidance issued by the RPSGB and maintain an internal temperature of between +2°C and +8°C. They are also equipped with a thermometer that records both the maximum and minimum temperatures.

Control and Efficiency

The neat, integral digital control panel on the outside of the fridge is easy to use and makes checking the internal temperature straightforward. All the pharmacy fridges have fans delivering optimum temperature stability and a faster temperature recovery after door openings. Automatic defrost ensures more efficient cooling.

Security

All Lec Medical pharmacy fridges come complete with door locks and two keys ensuring that only authorised individuals are able to access them. A temperature alarm, which sounds if the temperature deviates from the set parameters for more than 15 minutes, comes as standard across all pharmacy refrigerators. Open door sensors also come as standard across the pharmacy range. The auxiliary output* gives the user the option of connecting support systems such as a mains failure alarm signal.

Flexibility

Reversible doors* give greater flexibility where space is at a premium or where positioning the product is specific. Adjustable feet are also important when positioning the product on uneven surfaces. Each fridge comes with removable shelves giving maximum storage capacity.

Environmentally Friendly

Low energy consumption. CFC and HFC free.

Options

Lec Medical also offers several bespoke options across many pharmacy products, namely chart recorders, mains failure alarms, castors and European plugs. These are subject to an additional lead time and cost. For further details please contact our sales office on **0844 815 3742**.

PE102
SKU **444447386**
Product dimensions (mm) **379w x 500d x 702h**
Packed dimensions (mm) **450w x 558d x 740h**
Product weight (kg) **23.5**
Packed weight (kg) **25.5**
Net capacity (litres) **35**
Gross capacity (litres) **41**

PE207
SKU **444448564**
Product dimensions (mm) **502w x 540d x 660h**
Packed dimensions (mm) **520w x 595d x 675h**
Product weight (kg) **27.5**
Packed weight (kg) **30.5**
Net capacity (litres) **67**
Gross capacity (litres) **82**

PE507
SKU **444448508**
Product dimensions (mm) **558w x 600d x 863h**
Packed dimensions (mm) **600w x 610d x 900h**
Product weight (kg) **34**
Packed weight (kg) **35.5**
Net capacity (litres) **137**
Gross capacity (litres) **153**

PE907
SKU **444448510**
Product dimensions (mm) **595w x 600d x 1550h**
Packed dimensions (mm) **625w x 700d x 1610h**
Product weight (kg) **81.5**
Packed weight (kg) **84**
Net capacity (litres) **275**
Gross capacity (litres) **286**

PE357
SKU **444448512**
Product dimensions (mm) **595w x 600d x 1850h**
Packed dimensions (mm) **625w x 700d x 1909h**
Product weight (kg) **85**
Packed weight (kg) **88**
Net capacity (litres) **340**
Gross capacity (litres) **354**

PE1607
SKU **444448567**
Product dimensions (mm) **685w x 910d x 1688h**
Packed dimensions (mm) **770w x 950d x 1810h**
Product weight (kg) **130**
Packed weight (kg) **140**
Net capacity (litres) **393**
Gross capacity (litres) **444**

PE102



PE207



PE507



PE907



PE357



PE1607



*Not available on all products. Please refer to specifications on pages 36-39.