



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax:01759 301143

sales@wolflabs.co.uk



If the literature below doesn't answer all your questions, click these links to:

[Request details on an alternative product of the same type](#)

[Request details on accessories](#)

[Request details on other product categories](#)



PRODUCT INFORMATION

Model: RLDF0505

Type: Labcold Pharmacy Refrigerator

A purpose-built refrigerator for laboratory, surgery and general use for the storage of drugs and vaccines. Conforms in full with NEW UK Department of Health “Green Book” requirements, plus recommendations in Royal Pharmaceutical Society of Great Britain’s “Code of Ethics and Standards”.



Capacity: 150 litres

Construction: White enamel exterior, fully moulded plastic interior.

Shelves supplied: 3, plastic coated wire mesh.

Exterior dimensions (mm): 820H x 600W x 600D

Features:

- Integral digital temperature display**
- Max/min temperature record with reset**
- Audible alarm with mute facility**
- Fan assisted air circulation**
- Internal light**
- Door lock, with two keys**
- Integral door handle**
- CFC, HCFC and Ammonia free**
- Accurate temperature control: +2°C to +8°C**
- Automatic defrost, with condensate removal**

Order Reference: RLDF 0505

Optional Accessories: Pharmacy Drawer System,

Order Reference: PHARMDRAWER01