



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax:01759 301143

sales@wolflabs.co.uk



If the literature below doesn't answer all your questions, click these links to:

[Request details on an alternative product of the same type](#)

[Request details on accessories](#)

[Request details on other product categories](#)



PRODUCT INFORMATION

Model: RLDF1005, RLDG1005 (Glass Door)

Type: Labcold Pharmacy/Vaccine Refrigerator

A purpose-built refrigerator for laboratory, surgery and general use for the storage of drugs and vaccines. Conforms in full with NEW UK Department of Health “Green Book” requirements, plus recommendations in Royal Pharmaceutical Society of Great Britain’s “Code of Ethics and Standards”.



Capacity: 300 litres

Construction: White enamel exterior, plastic interior

Shelves Supplied: 6 wire mesh

Exterior Dimensions (mm): 1487H x 600W x 650D

Features:

Digital display of temperature with max/min record

Audible alarm with mute facility

Max/min temperature memory and display

Fan assisted air circulation

Internal light

Door lock

CFC and HFC free and Ammonia free

Accurate temperature control: +2°C to +8°C

Integrated digital temperature display and controller

Automatic Defrost

Order Reference: RLDF1005, RLDG1005 Glass Door