



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax:01759 301143

sales@wolflabs.co.uk



If the literature below doesn't answer all your questions, click these links to:

[Request details on an alternative product of the same type](#)

[Request details on accessories](#)

[Request details on other product categories](#)



PRODUCT INFORMATION

Model: RLDF18041

Type: Labcold Pharmacy/Vaccine Refrigerator

For laboratory, surgery and general use for the storage of drugs and vaccines. Conforms in full with New UK Department of Health "Green Book" requirements.



Image for illustrative purposes only.

Capacity: 505 litres (gross)

Construction: White exterior and interior

Shelves supplied: 5, plastic coated wire mesh

Exterior dimensions (mm): 1695H x 770W x 730D

Features:

Digital display of temperature with max/min record

Audible alarm with mute facility

Fan assisted air circulation

Door lock

CFC, HCFC and Ammonia free

Accurate temperature control: +2°C to +8°C

Automatic defrost

Order reference: RLDF 18041