



**Wolflabs**

# Wolf Laboratories Limited

[www.wolflabs.co.uk](http://www.wolflabs.co.uk)

Tel: 01759 301142

Fax:01759 301143

[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)



**If the literature below doesn't answer all your questions, click these links to:**

[Request details on an alternative product of the same type](#)

[Request details on accessories](#)

[Request details on other product categories](#)



## Cleaning and Disinfection

60cm wide Freestanding/built-under units

### Performance specifications

Narrow-necked glassware: 39/E 329

Pipettes: 116/E 406

Test tubes: 1600/E 103-139

### Control system

MTN Multitronic Novo Lab controls with 4 programmes

Temperature and countdown indicators

### Features

Stainless-steel cabinet with 2 levels

Spray arms/Stainless steel

Docking system for upper baskets/mobile injector units

Quadruple filtration system

UP Qmax circulation pump 200 l/min

AP Drain pump

WES Water softener

AD pump for demin water

Supplied with 13A moulded plug

### Heating

Electric (2.5kW)

### Dispensing technology

1 door dispenser for powder detergent

2 door dispenser for liquid products (Duo Dos)

1 DOS 10/30 dispenser pump for liquid acidic agent

DOS module connection: DOS G 60 dispenser module for liquid detergent,

DOS G 10 dispenser module for liquid neutralising

### Cleaning technology

Fresh water system

Max. cleaning temperature 70°C

### Design

Frontloading unit with bottom-hinged door

Freestanding with /built-under unit without lid

Slanted fascia panel

AE Outer panelling stainless steel

H x W x D: 850 (820) x 600 x 600 mm

### Test certificates

VDE, VDE-EMV, DVGW

G 7804

