



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax:01759 301143

sales@wolflabs.co.uk



If the literature below doesn't answer all your questions, click these links to:

[Request details on an alternative product of the same type](#)

[Request details on accessories](#)

[Request details on other product categories](#)





UF1380

Solid Door Freezer



Features and Benefits

- Static Cooling
- Reversible Door
- Adjustable Feet
- Easy Clean Interior
- Digital Temperature Display
- Matching Chiller Page 109
- Adjustable Temperature

Technical Specifications

Temperature Range	-18/-22°C
Refrigerant	R134A
Type of Defrost	Manual
No of Shelves	6 Fixed Static
Exterior Finish	White
Interior Finish	White
Power Required	13 Amp Supply

External Dimensions	Height:	Width:	Depth:	Capacity:	Weight:
UF1380	1835mm	595mm	600mm	372ltr. (13cu.ft.)	90kg
Internal Dimensions					
UF1380	1541mm	448mm	445mm		