



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax:01759 301143

sales@wolflabs.co.uk



If the literature below doesn't answer all your questions, click these links to:

[Request details on an alternative product of the same type](#)

[Request details on accessories](#)

[Request details on other product categories](#)



WOLF LABORATORIES

Colenso House, Deans Lane, Pocklington, York YO42 2PX

Tel: 01759 301142 Fax: 01759 301143 Email: sales@wolflabs.co.uk Website: www.wolflabs.co.uk

Features and Benefits

- Auto Defrost
- Fan Assisted Cooling
- Lock Fitted
- Fully Automatic Operation
- Adjustable Shelving
- Reversible Door (Factory option)
- Adjustable Feet
- Temperature Display
- Matching Freezer Page 148
- Glass door available Page 19

Technical Specifications

| | |
|-------------------|---------------|
| Temperature Range | +1/+12°C |
| Refrigerant | R134A |
| Type of Defrost | Off Cycle |
| No of Shelves | 7 |
| Exterior Finish | White |
| Interior Finish | White |
| Power Required | 13 Amp Supply |

SD1380

Solid Door Refrigerator

TEFCOLD 



| External Dimensions | Height: | Width: | Depth: | Capacity: | Weight: |
|----------------------------|---------|--------|--------|--------------------|---------|
| SD1380 | 1825mm | 595mm | 640mm | 372ltr. (13cu.ft.) | 115kg |
| <i>Internal Dimensions</i> | | | | | |
| SD1380 | 1550mm | 515mm | 470mm | | |