



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax:01759 301143

sales@wolflabs.co.uk



If the literature below doesn't answer all your questions, click these links to:

[Request details on an alternative product of the same type](#)

[Request details on accessories](#)

[Request details on other product categories](#)



WOLF LABORATORIES

Tel: 01759 301142 Fax: 01759 301143 Email: sales@wolflabs.co.uk
Website: www.wolflabs.co.uk

Cryo Boxes



Cardboard boxes WxD 133x133mm

All boxes listed below have an external measure of 133x133 mm. Our boxes are offered in 3 different heights, 50, 75 and 100 mm.

The boxes are made of fibreboard/cardboard and are Water-repellent for longer durability.

The table below shows the article number by height.

Height	Article no. Box 50 mm	Article no. Box 75 mm	Article no. Box 100 mm
White	TE22005	TE22050	TE22100



Grid dividers

For the boxes above we offer 13 different grid dividers. Choose the number of cells and for which box height the divider should be used.

The 30 mm is for 50 mm boxes
The 40 mm is for 75 mm boxes
The 65 mm is for 100 mm boxes

Grid dividers, height	Article no. 30 mm for 50 mm box	Article no. 40 mm for 75 mm box	Article no. 65 mm for 100 mm box
5x5 cells for vials up to 25 mm O.D.	TE22111	-	-
7x7 cells for vials up to 16 mm O.D.	TE22106	TE22130	TE22140
8x8 cells for vials up to 14 mm O.D.	TE22108	TE22132	TE22142
9x9 cells for vials up to 13 mm O.D.	TE22110	TE22134	TE22144
10x10 cells for vials up to 12 mm O.D.	TE22120	TE22136	TE22146