



Wolflabs

Wolf Laboratories Limited

www.wolflabs.co.uk

Tel: 01759 301142

Fax:01759 301143

sales@wolflabs.co.uk



If the literature below doesn't answer all your questions, click these links to:

[Request details on an alternative product of the same type](#)

[Request details on accessories](#)

[Request details on other product categories](#)



PhotoDoc-It™ Imaging Systems


Compact system combines a digital color camera, hood, transilluminator and a color printer for low cost prints. The optional LCD monitor mounts to the front of the hood for preview and capture of gels.

High resolution **digital color camera** to capture gel images

Lightweight hood, with handles for portability, fits over transilluminators with filter sizes up to 21x26cm

System can be configured with choice of **UV transilluminator**

The PhotoDoc-It is designed to **replace instant Polaroid photodoc systems**




Includes New Viewport Window!

Ethidium Bromide filter located inside the hood

Viewport Window is UV blocking and allows visibility of gels

Quickly generate **4 x 6" color prints** (108 prints)



- Illuminate stained gels with UV fluorescence of benchtop transilluminator
- Capture brilliant color images with the megapixel resolution digital camera
- Save images to the memory card for transfer to a computer (USB connection) for viewing or documentation
- Instantly print images on the lab quality color photo printer (included)

PhotoDoc-It Systems include: Digi Camera 55 (camera specifications subject to change without notice), Hood Enclosure, EtBr Filter, Combination Lock, Color Printer and Choice of Benchtop UV Transilluminator

Ordering Information

PhotoDoc-It 55 with color printer and M-20 Transilluminator 97-0274-04 (100-115V)
97-0274-05 (230V)

PhotoDoc-It 55 with color printer and LM-20 Transilluminator 97-0274-07 (100-115V)
97-0274-08 (230V)

PhotoDoc-It 55 with color printer and M-26 Transilluminator 97-0274-10 (100-115V)
97-0274-11 (230V)

PhotoDoc-It 55 with color printer and LM-26 Transilluminator 97-0274-13 (100-115V)
97-0274-14 (230V)

LCD Screen (optional, not shown)
6.4" TFT Video Monitor; 98-0093-01
12 VDC Power

PhotoDoc-It Specifications

Digi Camera 55: Megapixel resolution color camera* with USB connectivity, 4X optical zoom lens

Hood: Compact design of metal construction
Handles for portability
Fits up to 21x26cm filters

Transilluminator (choose from):
M-20: 302nm UV, 20x20cm filter
LM-20: 302/365nm UV, 20x20cm filter
M-26: 302nm UV, 21x26cm filter
LM-26: 302/365nm UV, 21x26cm filter

Filter: 50mm Sq. Ethidium Bromide

Printer: Color, with 300x300 resolution 4x6" prints (108 prints/ink included) and USB connectivity
Dim: 7 x 5 x 2.5 in. (178 x 127 x 64mm)

Combination Lock: 3 digit resettable combination lock

THERMAL PRINTER

Document gels, films and membranes with a high quality thermal printer. Easily connects for use with UVP's imaging systems (except PhotoDoc-It). Process low cost prints in seconds. Replacement paper rolls install easily.



Digital Thermal Printer for archive quality prints in seconds.

Dimensions:
10.74H x 13L x 3W in.
(273 x 330 x 76mm)

- Compact printer size takes up minimal lab space
- Supplies 256-gray scale for high resolution, archive-quality thermal film prints

Ordering Information

Digital Thermal Printer Mitsubishi 89-0069-06 (115V)
89-0069-07 (230V)

Digital Thermal Printer Sony (for Europe) 89-0357-01 (230V)

Thermal Paper, Mitsubishi 4 rolls (800 prints) 89-0038-01

Thermal Paper, Sony super glossy UPP110HG 5 rolls (1000 prints) 89-0174-01

Thermal Paper, Sony glossy UPP110HA 5 Rolls (1000 prints) 89-0031-01