

# WOLF LABORATORIES

Tel: 01759 301142 Fax: 01759 301143 Email: sales@wolflabs.co.uk  
Website: www.wolflabs.co.uk

## Memory & Alarm Thermo-hygrometers

Consistent & accurate temperature and humidity monitoring with 2 year warranty



These tried-and-tested thermo-hygrometers monitor the atmosphere every 18 seconds to ensure safe, dry conditions. With a world-renowned Testo™ humidity sensor built into the water-resistant housing, they will provide accurate measurements for years.

- Displays temperature in °C or °F
- Displays humidity in %RH
  - Large, clear screen
  - Dew point calculation
- Min/max value display with reset
- Monitors upper & lower limits with visual alarm (**AHT-608H2**)
  - Integrated patented Testo™ humidity sensor
  - Hang from wall or stand upright on a table
    - Low battery indicator
    - 2 year warranty
- **Supplied with battery**
- **Optional Calibration Certificate available**

ACAL-HT Calibration Certificate - Issued at 2 Temperature points: 20°C & 40°C and 2 RH Points: 35% & 75% @20°C. (Additional points are available if required. These must be specified at the time of order. Please call for price and details.)

### Technical Specifications

	<b>AHT-608H1</b>	<b>AHT-608H2</b>
Temperature Range	0 to 50°C / 32 to 122°F	-10 to 70°C / 32 to 122°F
Humidity Range	+10 to +95% RH	+2 to +98% RH
Dewpoint Parameters	-20 to 50°C	-40 to 70°C
Resolution	0.1°C/F (0 to +50°C/32 to 122°F)	0.1°C/F (-10 to +70°C/14 to 160°F)
Accuracy (± 1 digit)	±0.5°C/1°F at 25°C/77°F ±3% RH(+10 to +95% RH)	±0.5°C/1°F at 25°C/77°F ±2% RH(+2 to +96% RH)
Power	9V battery	
Dimensions	120 x 89 x 40mm	
Weight	168g	

### Related Product Codes

- ACAL-HT Humidity Calibration Certificate