


WOLF LABORATORIES

Tel: 01759 301142 Fax: 01759 301143 Email: sales@wolflabs.co.uk
Website: www.wolflabs.co.uk

Pipette Fillers



Quick-Release Pipette Fillers

Our pipetting aids feature a rapid release valve for one-handed operation. Stabilizing fingers in the collar of the product hold pipettes securely. They work with glass or plastic pipettes. You draw or release liquid by rotating the thumb wheel. To empty rapidly you merely apply light pressure on the fast release lever.

Cat No.	Description	Units
ZXZHS20619A	For pipettes up to 2mL	Each
ZXZHS20619B	For pipettes up to 10mL	Each
ZXZHS20619C	For pipettes up to 25mL	Each