

Med-Air

Article:1153000014

- Oil free compressor technology
- Maintenance free except for the intake filter(s) and sinter filter
- Ready for connection
- Contained within a temperature controlled acoustic cabinet with lockable doors
- Low noise level (lower than first generation units)
- Suitable for use in an S1 condition up to 30°C ambient
- Medical grade castors for ease of movement
- Auxiliary connection to enable external pressure vessel to be added if required



Technical data

Delivery Rate @ 5 bar	124 l/min
Receiver volume	12.5 liters
Voltage	230 V
Frequency	50/60 Hz
Current	8.6A
Protection class	IP44
Dimensions in mm (H x W x D)	945 x 670 x 455
Weight	102 Kg (Unpacked)



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.