



Product Code: XR1300H

xtra by Foster 1300 Litre Upright Refrigerated Cabinet




Features & Benefits

- The xtra range: Built for you, designed for your budget
- Reliable and robust refrigeration at an affordable price
- In stock and available now for immediate shipping
- Robust stainless steel finish
- Easy read microprocessor controller provides simple and accurate temperature control
- Full Foster 2 years parts and labour warranty

Shelving		Doors		Refrigeration	
6 shelves (GN2/1)	Standard	Left and right hinged door	Standard	R290 - Energy Efficient	Standard
Electrical & Service Data		Solid	Standard	R134a - Industry Standard	Option
230/50/1	Standard	Specification		Integral Refrigeration System	Standard
220/60/1 Electrical Supply	Option	Refrigerator +2/+8°C	Standard	Finish	
Fuse Required (A)	10	Width (mm)	1390	s/s ext & alu int	Standard
Refrigerant Charge (kg)	0.420	Depth (mm)	850		
Energy		Height (mm)	1985		
kWh (24 hours)	4.4565	Gross Volume (Ltr)	1300		
AEC	2027.94				
Nett Volume	901.2				
Energy Grade	G				
Global Warming Potential	1430				
Energy Efficiency Index	97.04				

Footnotes

- Extra shelf and pair of trayslides available as an option
- Supplied with side integrated door handles in lieu of recessed black plastic handles.



Application

This appliance is intended for use in ambient temperatures up to {1}°C

Safety

Precautionary instructions and maintenance advice, to ensure efficient operation, can be found in this products operating instructions.

Recycling

For recycling and end of life care for this product, please refer to environmental section of the website or visit [our policies](#)



Energy



Quality



Value



Design



Hygiene



After Care

All images shown are for illustration purpose only. All images and photographs are not intended to be relied upon for, nor to form part of, any contract unless specifically incorporated in writing into the contract. Solely the dimensions and specifications incorporated in the quotation or purchase order confirmation shall be contractually binding.



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.