

Haier

UFR106 Infrared Thermometer Introduction

额温枪
国内定制专用校准文件

UFR106-- Product Features

Haier



Measure
In 1s

Automatic
power-off

20
memories

Switch
between
Celsius and
Fahrenheit



Measure in 1 second

No contact with human skin, one second temperature measurement



Green Backlight

In low light, it can also clearly display the screen information, so that you can easily read the data.

UFR106--Product Parameters

Haier



ITEM	UFR106
Measurement method	Non-contact
Measurement range	32.0℃~43.0℃ (89.6℉~109.4℉)
Accuracy	±0.2℃/±0.4℉ (within 36℃~39℃ /96.8℉~102.2℉)
Display resolution	0.1℃/0.1℉
Memory function	20 sets memory of measurement values
Power source	2pcs AAA alkaline battery Notice: no battery
Main unit weight	Approx. 75g (batteries not included)
Main unit size	L151mm * W53mm * H41mm



ITEM	UFR106
BOX SIZE	93*53*176mm
CTN SIZE	392*290*376mm
CTN	40pcs
G.W	6.5KG
N.W	5.9KG





WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.