



LSF151UK ATEX Under-counter Laboratory Freezer Solid Door 151L

LSF151UK ATEX

- ✓ -18°C to -25°C temperature range
- ✓ External digital controller and temperature display
- ✓ Fitted lock and 2 keys
- ✓ Spark free interior
- ✓ Audible and visual temperature alarms

SPECIFICATION

Dimensions	595mm W x 845mm H x 595mm D
Weight	36kg
Product name	Under-counter Laboratory Freezer Solid Door 151L LSF151UK ATEX
GDPA code	LSF151UK ATEX
SKU code	UK 444442773, EU 444442774
External dimensions	595mm W x 845mm H x 595mm D
Internal dimensions	500mm W x 705mm H x 450mm D (step 185mm H x 165mm D)
Packed dimensions	657mm W x 872mm H x 694mm D
Product weight	36kg
Packed weight	40kg
Gross capacity	151L
Net capacity	151L
External finish	White
Colour	White/green
Internal finish	Plastic
Type of plastic	SRL800 (HIPS/ABS mix)
Door type	Solid
Door style	Hinged, reversible
Door open depth	1165mm
Replacement door seal	N
Open door sensor	Y
Lock	Y
Adjustable feet	1 adjustable front foot
Temperature range	-18°C to -25°C
Temperature display on panel	Y

Temperature alarm	Y
Number of baskets	3
Basket material	Wire
Number of shelf positions	Fixed
Interior light	N
Auto or manual defrost	Manual
Fan	N
Controller type	Digital, external
Climate class	SN
Refrigerant type	R600a
CFC free	Y
HFC free	Y
HCFC free	Y
UK 3 pin plug	Y
Voltage supply	230
Spark free	Y
Auxiliary output	Y
Insulation thickness	45mm
Insulation type	Cyclopentane
Air tempertaure alarm	Y
Door open alarm	Y
Stacking kit compatible	Y



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.