



## Freestanding Laboratory Refrigerator Solid Door 288L LSR288UK ATEX

LSR288UK ATEX

- ✓ 2°C to 10°C temperature range
- ✓ External digital controller and temperature display
- ✓ Fitted lock and 2 keys
- ✓ Spark free interior
- ✓ Audible and visual temperature alarms

## SPECIFICATION

|                     |                                                                       |
|---------------------|-----------------------------------------------------------------------|
| Dimensions          | 595mm W x 1570mm H x 634mm D                                          |
| Weight              | 58kg                                                                  |
| Product name        | Freestanding Laboratory Refrigerator Solid Door 288L<br>LSR288UK ATEX |
| GDPA code           | LSR288UK ATEX                                                         |
| SKU code            | UK 444440966 EU 444440967                                             |
| External dimensions | 595mm W x 1570mm H x 634mm D                                          |
| Internal dimensions | 480mm W x 1382mm H x 445mm D (step 255mm H x 175mm D)                 |
| Packed dimensions   | 657mm W x 1585mm H x 734mm D                                          |
| Product weight      | 58kg                                                                  |
| Packed weight       | 63kg                                                                  |
| Gross capacity      | 284L                                                                  |
| Net capacity        | 284L                                                                  |
| External finish     | White                                                                 |

|                              |                           |
|------------------------------|---------------------------|
| Colour                       | White/green               |
| Internal finish              | Plastic                   |
| Type of plastic              | SRL800 (HIPS/ABS mix)     |
| Door type                    | Solid                     |
| Door style                   | Hinged, reversible        |
| Door open depth              | 1185mm                    |
| Replacement door seal        | N                         |
| Open door sensor             | Y                         |
| Lock                         | Y                         |
| Adjustable feet              | 2 front, 2 rear rollers   |
| Temperature range            | 2°C to 10°C               |
| Temperature display on panel | Y                         |
| Temperature alarm            | Y                         |
| Number of shelves            | 5                         |
| Shelf material               | Wire                      |
| Number of shelf positions    | 12                        |
| Shelf size                   | 540mm W x 6mm H x 497mm D |
| Moulded drawer               | 1 solid                   |
| Adjustable door bins         | 4 solid                   |
| Interior light               | N                         |
| Auto or manual defrost       | Auto                      |
| Fan                          | N                         |
| Controller type              | Digital, external         |
| Climate class                | SN                        |
| Refrigerant type             | R600a                     |
| CFC free                     | Y                         |
| HFC free                     | Y                         |
| HCFC free                    | Y                         |
| UK 3 pin plug                | Y                         |

|                         |              |
|-------------------------|--------------|
| Voltage supply          | 230          |
| Spark free              | Y            |
| Auxiliary output        | N            |
| Insulation thickness    | 60mm         |
| Insulation type         | Cyclopentane |
| Air tempertaure alarm   | Y            |
| Door open alarm         | Y            |
| Stacking kit compatible | N            |



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.