



Countertop Laboratory Freezer Solid Door 50L ISU27C

ISU27C

- ✓ -18°C to -25°C temperature range
- ✓ External digital controller and temperature display
- ✓ Fitted lock and 2 keys
- ✓ Spark free interior
- ✓ Audible and visual temperature alarms

SPECIFICATION

Dimensions	500mm W x 562mm H x 510mm D
Weight	17kg
Product name	Countertop Laboratory Freezer Solid Door 50L ISU27C
GDPA code	ISU27C
SKU code	UK 444441802, EU 444442351, AU 444442352
External dimensions	500mm W x 562mm H x 510mm D
Internal dimensions	345mm W x 365mm H x 330mm D
Packed dimensions	520mm W x 590mm H x 595mm D
Product weight	17kg
Packed weight	20kg
Gross capacity	50L
Net capacity	38L
External finish	White
Colour	White/green

Internal finish	Plastic
Type of plastic	SRL800 (HIPS/ABS mix)
Door type	Solid
Door style	Hinged, reversible
Door open depth	976mm
Replacement door seal	N
Open door sensor	N
Lock	Y
Adjustable feet	N
Temperature range	-18°C to -25°C
Temperature display on panel	Y
Temperature alarm	Y
Number of shelves	1
Shelf material	Wire
Number of shelf positions	Fixed
Interior light	N
Auto or manual defrost	Manual
Fan	N
Controller type	Digital, external
Climate class	N
Refrigerant type	R600a
CFC free	Y
HFC free	N
HCFC free	Y
UK 3 pin plug	Y
Voltage supply	230
Spark free	Y
Auxiliary output	N
Insulation thickness	30mm

Insulation type	Cyclopentane
Noise level	<45dB
Energy consumption	0.58 kW/24h
Air temperature alarm	Y
Stacking kit compatible	N



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.