



Under-counter Laboratory Freezer Solid Door 102L

ISU37C

- ✓ -18°C to -25°C temperature range
- ✓ External digital controller and temperature display
- ✓ Fitted lock and 2 keys
- ✓ Spark free interior
- ✓ Audible and visual temperature alarms

SPECIFICATION

Dimensions	495mm W x 830mm H x 620mm D
Weight	33kg
Product name	Under-counter Laboratory Freezer Solid Door 102L ISU37C
GDPA code	ISU37C
SKU code	UK 444441803, EU 444442353, AU 444442354
External dimensions	495mm W x 830mm H x 620mm D
Internal dimensions	414mm W x 683mm H x 415mm D
Packed dimensions	560mm W x 850mm H x 655mm D
Product weight	33kg
Packed weight	36.5kg
Gross capacity	102L
Net capacity	87L
External finish	White
Colour	White/green
Internal finish	Plastic
Type of plastic	SRL800 (HIPS/ABS mix)
Door type	Solid
Door style	Hinged, reversible
Door open depth	1050mm
Replacement door seal	N
Open door sensor	N
Lock	Y
Adjustable feet	1 adjustable front foot
Temperature range	-18°C to -25°C
Temperature display on panel	Y
Temperature alarm	Y

Number of shelves	2
Shelf material	Wire
Number of shelf positions	Fixed
Interior light	N
Auto or manual defrost	Manual
Fan	N
Controller type	Digital, external
Climate class	N
Refrigerant type	R134a
CFC free	Y
HFC free	N
HCFC free	Y
UK 3 pin plug	Y
Voltage supply	230
Spark free	Y
Auxiliary output	N
Insulation thickness	35mm
Insulation type	Cyclopentane
Noise level	<43dB
Energy consumption	1.3 kW/24h
Air temperature alarm	Y
Stacking kit compatible	N



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.