



British healthcare and laboratory
refrigeration specialists



Laboratory Range

Under Counter Fridge and Freezer



FEATURES	KEY BENEFITS
Spark free interior	Our range of laboratory refrigerators and freezers offers the safe and reliable storage of volatile liquids, reagents, chemicals, poisons and samples with a spark-free interior that reduces the risk of internal explosion.
Lock	Security is essential in laboratories so our Lec Medical fridge and freezer come complete with door lock and two keys, ensuring that only authorised individuals are able to access them.
12 month warranty	Our under counter fridge and freezer come complete with a 12 month warranty for added peace of mind.




British healthcare and laboratory
refrigeration specialists

Laboratory Range

Under Counter Fridge and Freezer



SPECIFICATION TABLE



Internal Drawer



Lock

Product Description

Product Dimensions (external)

mm

Capacity

litres

Colour/Finish

White

Reversible Door

Y/N

Yes

Door Lock

Y/N

Yes

Drawer Number

1 Plastic

Shelf Number

3 Glass

Handle Type

Recessed

Defrost Capability

Auto Defrost

Annual Energy Consumption

kwh

117

Thermostat

Adjustable

Spark Free

Y/N

Yes

HFC Free

Y/N

Yes

CFC Free

Y/N

Yes

Product Description

Product Dimensions (external)

mm

Capacity

litres

Colour/Finish

White

Reversible Door

Y/N

Yes

Door Lock

Y/N

Yes

Drawer Number

3 Plastic

Shelf Number

1 Plastic

Handle Type

Recessed

Annual Energy Consumption

kwh

157

Thermostat

Adjustable

Spark Free

Y/N

Yes

HFC Free

Y/N

Yes

CFC Free

Y/N

Yes

Laboratory Fridge

545 x 554 x 845mm (W x D x H)

130

Laboratory Freezer

545 x 554 x 845mm (W x D x H)

85



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.